

mAb Charge Variant Analysis by Weak Cation Exchange (WCX) Chromatography



Monoclonal antibodies (mAbs) are a prominent biotherapeutic that are produced by most major pharmaceutical companies. mAbs are about 150,000 Daltons and have several modifications that add to the complexity of this class of macromolecule. Modifications can occur as the protein is produced or as they are manufactured and stored. Modifications such as C-terminal truncation, deamidation, changes in glycosylation, and amino acid deletion or substitution affect the overall charge of a mAb and are considered charge variants. Charge variants are typically monitored via ion exchange chromatography.^{1,2} Monitoring critical quality attributes (CQAs) such as charge variants is essential to ensure efficacy, safety, and immunogenicity of the final product is not compromised. Government agencies such as the U.S. Food and Drug Administration (FDA) require charge variant data when submitting a Biologics License Application.

Although charge variant analysis is an established technique, challenges remain for this assay. The buffers that are frequently used have a high salt concentration and can be corrosive to traditional stainless-steel hardware in liquid chromatography systems and columns. Using a bio-compatible or bio-inert system like the Agilent 1290 Infinity II Bio LC or the Agilent 1260 Infinity II Bio-Inert system mitigates this issue. The flow path is completely iron and stainless-steel free. All fittings and capillaries are made of an alloy, MP35N, so corrosion is reduced. In addition, this prevents reactions from corrosion such as oxidation from occurring. The 1290 Infinity II Bio LC is available as a binary pump which is ideal for producing accurate and precise gradients, especially for shallow gradients which are typical of charge variant analyses. The quaternary pump option is also useful for charge variant analysis for testing different buffer conditions which is made easier with Agilent Buffer Advisor software.

The Agilent Bio MAb column is a weak cation exchange column specifically designed to characterize the charge heterogeneity of monoclonal antibodies.³ It is available in PEEK hardware, which maintains an iron-free flow path throughout the 1290 Infinity II Bio LC system, essential for highly reproducible results (Figure 1).⁴

Furthermore, the column itself was specifically designed for mAb separations and has an extra dense and highly uniform weak cation exchange layer that is bonded to the hydrophilic polymeric coating, which is designed to eliminate non-specific interactions, making it ideal for biosimilar comparability studies (Figure 2).⁵




Figure 1. Acidic and basic variants of trastuzumab. Seven subsequent runs yields highly reproducible results in retention time and peak area.




Figure 2. Comparison of charge variants between rituximab and two biosimilars.

Best Practices for Optimizing Chromatographic Conditions

Two parameters are essential to achieve the desired resolution and enable optimal separation of charge variants: determining the optimal pH of the mobile phase as well as the optimal gradient slope. Both factors can have a major impact on the separation. Table 1 suggests parameters from where to start optimization.

	Salt Gradient	pH Gradient
Parameter	Value	Value
Column	Bio MAb, NP5, 2.1 x 250 mm, PEEK (part number 5190-2411)	Bio MAb, NP5, 2.1 x 250 mm, PEEK (part number 5190-2411)
Suggested LC System	Agilent 1290 Infinity II Bio LC System with High Speed Pump	Agilent 1290 Infinity II Bio LC System with Flex (Quaternary) Pump
Mobile Phase	A: 30 mM phosphate buffer, pH 6.8 B: 30 mM phosphate buffer, pH 6.8, 500 mM NaCl	A: Water B: 1.6 M NaCl C: 100 mM NaH ₂ PO ₄ D: 100 mM Na ₂ HPO ₄
Gradient	0-100 mM NaCl from 0-30 minutes; 100-500 mM from 30-31 min; Isocratic 500 mM from 31-35 min. Post-time 15 minutes	pH 6.0 to 8.0, 0 to 20 minutes 0 to 800 mM NaCl, 20 to 25 minutes 800 mM NaCl, 25 to 30 minutes
Flow Rate	0.25 mL/min	0.25 mL/min
Injection Volume	1-5 µL	1-5 µL
Temperature	Ambient	Ambient
Detection	280 nm	280 nm

Table 1. Suggested starting conditions

Mobile Phase Considerations

- The mobile phase should contain buffer to maintain the desired operating pH as well as consistent charge, the concentration is typically 20-30 mM.
- Phosphate buffers are very common in the range of pH 6 to 7. Other compatible buffers include acetate, Tris, and MES-containing buffers as well as acetonitrile and methanol. mAb samples must be soluble in the mobile phase.
- A competing ion must be introduced to elute the mAb from the column which is typically accomplished with 100 mM to 500 mM sodium chloride.
- Addition of sodium chloride will alter the pH of the mobile phase, so the pH will need to be readjusted.
- For optimal results, buffers should be made fresh and stored in the refrigerator unless being actively used as bacterial growth is common in dilute buffer. In addition, it is best practice to filter buffers to prevent column clogging.

pH Gradient Buffer Considerations

- Salt gradients are more common in ion-exchange chromatography; however, pH gradients are an alternate method^{6,7} that may result in higher resolution. The Bio MAb column is stable from pH 2-12.
- mAb pH gradients typically start at pH 6 and run up to pH 7-8 and have a sodium chloride salt column cleanup/ equilibration step.
- Monosodium phosphate and disodium phosphate are typical elution buffers.

More General Protein Mobile Phase Considerations

- The pH of the starting buffer should be 0.5 to 1 pH unit below the proteins' pI for cation-exchange chromatography.
- pH 6 is a good starting place if the pI of the protein is not known.

Column Dimensions

- The inner diameter (id) of the column should be selected based on the amount of mAb being analyzed. A general guideline is that the volume injected should not exceed more than 1-2% of the entire column volume. For example, a 2.1 x 50 mm column should ideally have 1.7 μ L or less injected while a 4.6 x 50 mm could handle an 8.3 μ L injection. If a 2.1 mm id column is selected, the LC should be plumbed to minimize dispersion.

- A longer column such as 250 mm will yield higher resolution, while a shorter column like 50 mm can help reduce run time.

Particle Size

- Generally, smaller particle sizes provide more efficient separation, but have higher operating pressures.
- Larger biomolecules like mAbs have a slower rate of diffusion, so particle sizes do not have as strong of an effect on resolution.
- Eluents with aqueous buffers are relatively viscous and can yield higher back pressures.
- Bio MAb particle sizes and max back pressures are as follows:
 - 1.7 μ m: 689 bar
 - 3 μ m: 551 bar
 - 5 μ m: 413 bar
 - 10 μ m: 275 bar
- Another consideration is that PEEK hardware has a pressure limit of 400 bar and is only available in the 5 and 10 μ m particle sizes. Salt gradients can be corrosive to stainless steel, so using PEEK columns in addition to a bio-inert or bio-compatible LC system will maintain robustness of the system.

Column Lifetime and Reproducibility

- Flow rates range from 0.1- 1.0 mL/min. 2.1 mm id columns typically run at 0.2-0.4 mL/min while 4.6 mm id columns typically run from 0.5-1.0 mL/min. Backpressure should be monitored and initially flow rates should be ramped slowly to ensure the column does not exceed the operating pressure limit.
- Bio MAb columns can withstand up to 80°C. However, to extend column lifetime, they should be regularly operated between 10-50°C.
- For maximum reproducibility, the equilibration/cleanup step in the gradient should be 5-10 column volumes.
- Consider using a guard column to extend the life of the analytical column.
- For extended storage, the column should be flushed for at least 15 column volumes and stored in 20 mM phosphate buffer with 0.1% sodium azide at pH 6.

Easy selection and ordering information

To order items listed in the tables below from the Agilent online store, add items to your Favorite Products list by clicking on the MyList # header links. Then, enter the quantities for the products you need, Add to Cart and proceed to checkout. Your list will remain under Favorite Products for your use with future orders.

If this is your first time using Favorite Products, you will be asked to enter your email address for account verification. If you have an existing Agilent account, you will be able to log in. However, if you don't have a registered Agilent account, you will need to register for one. This feature is valid only in regions that are e-commerce enabled. All items can also be ordered through your regular sales and distributor channels.

Description	Part No.
MyList 1: Bio MAb columns for charge variant analysis	
Agilent Bio MAb, NP5, 2.1 x 50 mm, PEEK	5190-2412
Agilent Bio MAb, NP5, 2.1 x 250 mm, PEEK	5190-2411
Agilent Bio MAb, NP5, 4.6 x 50 mm, PEEK	5190-2408
Agilent Bio MAb, NP5, 4.6 x 250 mm, PEEK	5190-2407
Agilent Bio MAb, NP10, 2.1 x 50 mm, PEEK	5190-2420
Agilent Bio MAb, NP10, 2.1 x 250 mm, PEEK	5190-2419
Agilent Bio MAb, NP10, 4.6 x 50 mm, PEEK	5190-2416
Agilent Bio MAb, NP10, 4.6 x 250 mm, PEEK	5190-2415
Agilent Bio MAb, NP1.7, 4 x 10 mm, guard	5190-2402
Agilent Bio MAb, NP3, 4 x 10 mm, guard	5190-2404
Agilent Bio MAb, NP5, 4 x 10 mm, guard	5190-2406
Agilent Bio MAb, NP10, 4 x 10 mm, guard	5190-2414
Agilent Bio MAb, NP1.7, 4.6 x 50 mm	5190-2401
Agilent Bio MAb, NP3, 4.6 x 50 mm	5190-2403
Agilent Bio MAb, NP5, 4.6 x 250 mm	5190-2405
Agilent Bio MAb, NP10, 4.6 x 250 mm	5190-2413
Agilent Bio MAb, NP5, 10 x 250 mm	5190-6884
Agilent Bio MAb, NP5, 21.2 x 250 mm	5190-6885
MyList 2: Standards	
Agilent NIST mAb, 25 µL	5191-5744
Agilent NIST mAb, 4 x 25 µL	5191-5745
MyList 3: Supplies & Solvents	
Connections & Tubing	
InfinityLab Quick Connect LC fitting	5067-5965
Agilent InfinityLab Quick Turn Fitting (for connection on column outlet)	5067-5966
Quick Turn Capillary MP35N 0.12 x 200 mm	5500-1595
Fitting Union Bio compatible, MP35N 10-32 coned. zero dead volume	5023-2625
Quick Turn Capillary MP35N 0.12 x 280 mm	5500-1596
Mounting tool for quick turn fittings	5043-0915
Inline pressure relief valve kit (For use when another detector is used in series after the fluorescence flow cell)	G4212-68001
6-Column Selector Valve, analytical, Bio-compatible. 1290 Infinity II Bio LC System	5320-0025
Quick Change 2-Position/10-Port Bio Valve Includes 1300 bar Quick Change Bio valve head.	5067-6682
Ultralow dispersion tubing kit for Agilent 1290 Infinity II Bio	5004-0007
Capillary kit for 2 position/10 port switching valve bio	5013-0002
Capillary Kit for 6-Column Selector Valve. Bio-compatible, 0.12 mm id 1290 Infinity II Bio LC System	5005-0070

Description	Part No.
MyList 3: Supplies & Solvents	
Sample Containment	
High recovery vial, screw top, with fixed insert, clear, 300 µL insert volume, 100/pk. Vial size: 12 x 32 mm (12 mm cap)	5188-6591
Cap, screw, blue, PTFE/red silicone septa, 100/pk. Cap size: 12 mm	5182-0717
Vial, crimp/snap top, polypropylene, 250 µL, 1,000/pk. Vial size: 12 x 32 mm (11 mm cap)*	5190-3155
Cap, snap, clear, PTFE/silicone/PTFE septa, 100/pk. Cap size: 11 mm (for 5190-3155)	5182-0566
InfinityLab 96-well plate, 0.5 mL, 30/pk	5043-9310
InfinityLab 96-well plate closing mat, 50/pk	5042-1389
Solvents & Additives	
InfinityLab Ultrapure LC/MS Water, 1 L	5191-4498
Solvent Filtration	
InfinityLab solvent filtration assembly	5191-6776
InfinityLab solvent filtration flask, glass, 2 L	5191-6781
Filter membrane, Nylon 47 mm, pore size 0.2 µm, 100/pk	5191-4341
Filter membrane, Regenerated Cellulose 47 mm, pore size 0.2 µm, 100/pk	5191-4340
Solvent bottle glass filter, solvent inlet, 20 µm	5041-2168
Solvent Handling	
InfinityLab solvent bottle, clear, 1 L	9301-6524
InfinityLab solvent bottle, amber, 1 L	9301-6526
Solvent bottle, clear, 2 L	9301-6342
Solvent bottle, amber, 2 L	9301-6341
InfinityLab waste can, GL45, 6 L with Stay Safe cap	5043-1221
InfinityLab charcoal filter with time strip, 58 g	5043-1193
Stay Safe starter kit and purging bottle, includes InfinityLab Stay Safe purging bottle (PN 5043-1339) and Stay Safe caps starter kit (PN 5043-1222)	5043-1340
InfinityLab Stay Safe cap starter kit	5043-1222
InfinityLab Stay Safe Purging Bottle	5043-1339

References:

1. Characterize mAb Charge Variants by Cation-Exchange Chromatography
[5991-5273EN](#)
2. Charge Variant Analysis, Application Compendium
[5994-2074](#)
3. Analysis of Intact and C-terminal Digested IgG1 on an Agilent Bio MAB 5 µm Column
[5991-0895EN](#)
4. How Shallow Can You Go? Refining charge variant analysis of mAbs with the Agilent 1290 Infinity II Bio LC System
[5994-2692EN](#)
5. Charge Variant and Aggregation Analysis of Innovator and Biosimilars of Rituximab
[5994-1496EN](#)
6. pH Gradient Elution for Improved Separation of Monoclonal Antibody Charge Variants
[5990-9629EN](#)
7. High-resolution Analysis of Charge Heterogeneity in Monoclonal Antibodies Using pH-gradient Cation Exchange Chromatography
[5991-1407EN](#)

Learn more:

www.agilent.com/chem/advancebio

Find a local Agilent customer center
in your country:

www.agilent.com/chem/contactus

U.S. and Canada

1-800-227-9770

agilent_inquiries@agilent.com

Europe

info_agilent@agilent.com

Asia Pacific

inquiry_lsca@agilent.com

DE99980273

This information is subject to change without notice.

© Agilent Technologies, Inc. 2023
Printed in the USA, May 10, 2023
5994-6053EN

