



crawford scientific

# GC Method Development


A Decision Tree for  
Gas Chromatography Method Development

---

BUILD BETTER GC METHODS WITH CRAWFORD SCIENTIFIC

# A DECISION TREE FOR GAS CHROMATOGRAPHY METHOD DEVELOPMENT


## Stage 1 – Is analyte volatile enough for GC analysis?



\*Remember that analytes often elute below their boiling point as the high linear gas flow reduces the vapour pressure within the system. Even if the inlet is at 250°C this may be enough to vaporise analytes with boiling points much higher than this.

# A DECISION TREE FOR GAS CHROMATOGRAPHY METHOD DEVELOPMENT

## Stage 2 – Nature of the sample and the requirement for sample preparation



\*Check the chemical compatibility of filter materials and housings when filtering organic solvents. Some assessment should also be made of the absorptivity of the filter material towards target analytes to ensure good accuracy and reproducibility.

# A DECISION TREE FOR GAS CHROMATOGRAPHY METHOD DEVELOPMENT

## Stage 3 (I) – Select appropriate detection technique: Quantitative Analysis

FID*	ECD	TCD	NPD	FPD	CLD
<b>Default choice</b>	<b>Analyte contains Oxygen / Halogens</b>	<b>Permanent Gas Analysis</b>	<b>Analyte contains N / P</b>	<b>Analyte contains S / P</b>	<b>Analyte contains S / N</b>
<ul style="list-style-type: none"> <li><b>Response Type:</b> Universal (C)</li> <li><b>Responds To:</b> Carbon</li> <li><b>Mass or Concentration Dependant:</b> Mass</li> <li><b>Destructive:</b> Y</li> <li><b>LOD:</b> <math>10^{-12}</math> gC/sec</li> <li><b>Linear Range:</b> <math>10^7</math></li> <li><b>Relative Selectivity:</b> none</li> <li><b>Uses:</b> Default detector choice for organic compounds</li> <li><b>Optimise:</b> <b>Temperature</b> (typically 300 – 350°C) <b>Air to hydrogen ratio</b> (typically 100:1) <b>Make-up gas flow</b> (match air flow and optimise in +/-20% steps)</li> </ul>	<ul style="list-style-type: none"> <li><b>Response Type:</b> Selective</li> <li><b>Responds To:</b> Halogen or Oxygen containing analytes</li> <li><b>Mass or Concentration Dependant:</b> Concentration</li> <li><b>Destructive:</b> N</li> <li><b>LOD:</b> <math>10^{-14}</math> g/mL</li> <li><b>Linear Range:</b> <math>10^4</math></li> <li><b>Relative Selectivity:</b> Up to <math>10^5</math> versus carbon</li> <li><b>Uses:</b> Environmental analysis for pesticides &amp; herbicides</li> <li><b>Optimise:</b> <b>Make up Gas Type &amp; Flow</b> (Nitrogen, Methane or Methane with 5% Argon) <b>Temperature</b> (250 - 350°C typical) <b>Current &amp; Mode</b> Magnitude of applied current / pulsed or constant current mode</li> </ul>	<ul style="list-style-type: none"> <li><b>Response Type:</b> Universal (C)</li> <li><b>Responds To:</b> Thermal Conductivity</li> <li><b>Mass or Concentration Dependant:</b> Concentration</li> <li><b>Destructive:</b> N</li> <li><b>LOD:</b> <math>10^{-9}</math> g/mL</li> <li><b>Linear Range:</b> up to <math>10^5</math></li> <li><b>Relative Selectivity:</b> none</li> <li><b>Uses:</b> Permanent gas analysis or where non-destructive general detector is required</li> <li><b>Optimise:</b> <b>Temperature</b> (250 - 350°C typical) <b>Bead Power / Voltage</b> (follow manufacturers guidelines) <b>Hydrogen Flow</b> (low enough not avoid a flame formation)</li> </ul>	<ul style="list-style-type: none"> <li><b>Response Type:</b> Universal (C)</li> <li><b>Responds To:</b> Nitrogen or Phosphorous</li> <li><b>Mass or Concentration Dependant:</b> Mass</li> <li><b>Destructive:</b> Y</li> <li><b>LOD:</b> <math>10^{12}</math> gN/sec</li> <li><b>Linear Range:</b> <math>10^5</math></li> <li><b>Relative Selectivity:</b> <math>2 \times 10^4</math> v C (N) <math>7 \times 10^4</math> v C (P)</li> <li><b>Uses:</b> Environmental analysis / Petrochemical analysis</li> <li><b>Optimise:</b> <b>Temperature</b> (lowest works best) <b>Reference flow</b> (typically 3x column + makeup flow) <b>Filament resistance</b> (follow manufacturer guidelines)</li> </ul>	<ul style="list-style-type: none"> <li><b>Response Type:</b> Universal (C)</li> <li><b>Responds To:</b> Phosphorous or Sulphur</li> <li><b>Mass or Concentration Dependant:</b> Mass</li> <li><b>Destructive:</b> Y</li> <li><b>LOD:</b> <math>10^{-13}</math> gP/sec <math>10^{-12}</math> gS/sec</li> <li><b>Linear Range:</b> <math>10^4</math> (P) Non-Linear (S)</li> <li><b>Relative Selectivity:</b> <math>10^6</math> v C (P) <math>10^6</math> v C (S)</li> <li><b>Uses:</b> Any analysis requiring trace determination of P or S containing compounds</li> <li><b>Optimise:</b> <b>Temperature</b> (up to 250°C typical) <b>Filter</b> (select correct filter- either P or S) <b>As per FID</b></li> </ul>	<ul style="list-style-type: none"> <li><b>Response Type:</b> Universal (C)</li> <li><b>Responds To:</b> Nitrogen or Sulphur</li> <li><b>Mass or Concentration Dependant:</b> Mass</li> <li><b>Destructive:</b> Y</li> <li><b>LOD:</b> <math>10^{-12}</math> gS/sec <math>10^{-12}</math> gN/sec</li> <li><b>Linear Range:</b> <math>&gt;10^4</math> (S &amp; N)</li> <li><b>Relative Selectivity:</b> <math>10^7</math> v C (S) <math>10^7</math> v C (N)</li> <li><b>Uses:</b> Petroleum characterisation, flavour analysis, environmental monitoring</li> <li><b>Optimise:</b> <b>Temperature Furnace Temperature</b> (800°C is typical) <b>Ozone Flow</b> <b>Hydrogen Flow</b> <b>Air Flow</b> (follow manufacturer guidelines)</li> </ul>

Move to next stage

\*FID – Flame Ionisation Detector / ECD – Electron Capture Detector / TCD – Thermal Conductivity Detector / NPD – Nitrogen Phosphorous Detector / FPD – Flame Photometric Detector / CLD – Chemiluminescence Detector. Also available – Photoionisation detector (PID) / Atomic Emission Detector (AED) / Electrolytic Conductivity Detector (ELCD) / Infrared Detector (IRD)

# A DECISION TREE FOR GAS CHROMATOGRAPHY METHOD DEVELOPMENT

## Stage 3 (II) – Select appropriate detection technique: Qualitative and / or Quantitative Analysis

### Mass Spectrometric Detection Electron Ionisation (EI) Mode

- **Response Type:**  
*Selective & Universal*
- **Responds To:**  
*Radical Cations created in ion source*
- **Mass or Concentration Dependant:**  
*Mass*
- **Destructive:**  
*Y*
- **LOD:**  
 *$10^6$*
- **Linear Range:**  
 *$10^7$*
- **Relative Selectivity:**  
*0.1Da in selected ion or higher in single or multiple reaction monitoring (SRM / MRM) modes*
- **Uses:**  
*Universal for qualitative analysis*
- **Optimise:**  
*Temperature of Transfer Line, Ion Source and Mass Analyser*  
*The MSD is a complex detector which requires both calibration and optimisation using onboard tuning algorithms*

Poor Molecular Ion Formation or  
Sensitivity in EI mode

Try EI at 25eV electron energy to reduce degree of fragmentation and promote intensity of molecular ion

Mass Spectrometric Detection  
Chemical Ionisation (CI) Mode


+ve      -ve

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>• <b>Chose Reagent Gas:</b><br/><i>Isobutane</i><br/><i>Methane</i><br/><i>Ammonia</i><br/><i>(in order of decreasing molecular ion intensity)</i></li> <li>• <b>LOD:</b><br/><i><math>10^{-6} g</math></i></li> </ul> | <ul style="list-style-type: none"> <li>• <b>Chose Reagent Gas:</b><br/><i>Isobutane</i><br/><i>Methane</i><br/><i>Ammonia</i><br/><i>(in order of decreasing molecular ion intensity)</i></li> <li>• <b>LOD:</b><br/><i><math>10^{-15} g</math></i></li> </ul> |
|---|--|

Move to next stage


# A DECISION TREE FOR GAS CHROMATOGRAPHY METHOD DEVELOPMENT

## Stage 4 – Select appropriate sample introduction technique<sup>1,2</sup>



# A DECISION TREE FOR GAS CHROMATOGRAPHY METHOD DEVELOPMENT

## Stage 5 – Select appropriate GC Column



# A DECISION TREE FOR GAS CHROMATOGRAPHY METHOD DEVELOPMENT

## Stage 6 – Develop & Optimise Oven Temperature Program <sup>1,2</sup>

### Splitless Injection

#### Set Initial Gradient Temperature & Hold Time

Follow screen method for Split injection but insert initial hold time of 1 min.

Maintain Initial Temperature at 20°C below boiling point of sample diluent

Optimise hold time empirically in steps of +/- 10 seconds until peak area of early eluting peaks is <1% RSD

Follow Split Temperature Program Optimisation

### Split Injection

Screen sample using the following conditions

Initial Temp.: 40°C  
 Initial Time: 0 mins.  
 Programming Rate: 10°C/min.  
 Final Temp.: 330°C  
 (or column gradient max. temp.)  
 Final Time: 10 mins

#### Peaks elute in a window of 5 mins or less?

Y

N

**Isothermal Analysis Possible**  
 Calculate  $T_{(f)}$  (retention temperature of final peak in the screening chromatogram)

$$T_{iso} = T_{(f)} - 45^{\circ}\text{C}$$

#### Set Initial Gradient Temperature

Calculate  $T_{(i)}$  (retention temperature of first analyte peak in the screening chromatogram)

$$T_{initial} = T_{(i)} - 45^{\circ}\text{C}$$

#### Optimise Temperature Ramp Rate

Calculate or measure  $t_0$  (hold-up time)  
 Optimum Ramp Rate = 10°C per  $t_0$   
 Further optimise in steps of +/- 10°C/min.

#### Optimise Final Temperature & Time

Final temperature is  
 $T_{(final)} = T_{(f)} + 20^{\circ}\text{C}$   
 Final Hold Time is  $t_0 \times 5$  mins (optimise empirically)

#### Optimise resolution for critical pairs using mid-ramp holds

Determine elution temperature of critical pair  $T_{(crit)}$  and insert mid ramp hold at  
 $T_{(iso-hold)} = T_{(crit)} - 45^{\circ}\text{C}$

Optimise the mid-ramp hold time empirically

1. The various relationships shown here are 'rules' of thumb' and need to be used as such to guide next steps. There are more finessed calculations which can be found at [www.CHROMacademy.com](http://www.CHROMacademy.com)
2. It is likely that some empirical optimisation will be required after each rule of thumb is employed and the results of the calculations will not deliver the optimum possible value for that variable.



# Crawford Scientific Ltd.

Holm Street, Strathaven  
ML10 6NB  
Scotland

## Contact Us

phone: 01357 522 961

email: [enquiries@crawfordscientific.com](mailto:enquiries@crawfordscientific.com)

web: [www.crawfordscientific.com](http://www.crawfordscientific.com)



## ABOUT THE AUTHOR

Tony Taylor has over 30 years of experience in developing chromatographic methods. As the Technical Director of three varied contract and application development laboratories, he understands what frustrates analytical chemists and how to help them overcome problems. He has helped thousands of budding chromatography method developers using his own experiences and insights, working with students to improve knowledge and understanding of chromatographic processes and their application.