


---

## 7694 Oven Insulation

	<b>Description</b>	<b>HP Part No.</b>	<b>Qty</b>
1	Oven, top cover	301-111-HSP	1
2	Thermalite	301-150-HSP	1
3	Durablanket th .38mm	301-151-HSP	1
4	Rod, lever lock cover	301-019-HSP	1
5	Carousel See figure 6-8		1
6	Frame, base case	301-110-HSP	1
7	Insulation cover	301-023-HSP	1
8	Disk, cover	301-120-HSP	1
<b>Associated screws</b>			
A	Screw, M4 X 16 flat head	0515-1014	4
B	Standoff, M5	300-914-HSP	2



**Figure 6-3.** Insulation