

## Flame Photometric Detector Covers and Electronics

Description	HP Part No.	Qty.
1. Flame photometric detector assembly		1
Single wavelength FPD	G1535-60550	
Dual wavelength FPD (not shown)	G1535-60560	
2. Detector top cover, single FPD (not shown)	G1535-80550	1
2A. Detector top cover, dual FPDs	G1535-80560	1
3. Pneumatics control module	G1535-60520	1
4. Electronics cover pan	G1535-00120	1
5. Electronics top cover	G1535-80540	1
Cover assembly (6890) (not shown)	G1535-80520	1
6. Actuator solenoid bracket (optional)	G1580-00070	–
7. Printed circuit board	G1535-60010	1




Figure 1320-14 Exploded Diagram of the FPD Covers and Electronics