

## Microcell Electron Capture Detector ( $\mu$ ECD) Manifold Assembly

Description	HP Part No.	Qty
1. EPC Control Valve (replacement kit)	G1531-60610	1
2. O-Rings HP M8	0905-1493	13
3. Screw, M4 x 8 mm, Torx, T20	0515-0910	1
4. EPC Micro ECD Manifold Assembly	G2397-60550	1
NS Screw, M3 x 25 mm, Torx, T10	0515-0683	8
NS Screw, M3 x .5 x 8 mm, Torx, T10	0515-0655	2
5. ECD Make-up Gas Assembly	G2397-20540	1
Chemical Traps	G1533-80540	2

Microcell Electron Capture Detector ( $\mu$ ECD) Manifold Assembly




Figure 1320-10 ( $\mu$ ECD) Manifold Assembly