

## Thermal Conductivity Detector (TCD) Manifold Assembly

Description	HP Part No.	Qty
1. EPC Control Valve (Order No. G1531-60610)	G1531-60541	2
1A. Non-EPC Detector On-Off Valve	G1531-60530	2
2. O-Rings HP M8 (Order No. 5180-4181, 12/pkg)	0905-1039	9
3. Screw, M4 x 6 mm, Torx, T20	0515-2832	2
4. EPC Control Module Assembly	G1532-60520	1
5. Non-EPC Control Module Assembly	G1532-60500	1
6. Screw, M3 x 8 mm, Torx T10	0515-0655	2
7. Make-up Gas, Manual Forward Pressure Regulator Assembly	G1532-60530	1
7A. Pressure Regulator, 0–100 psi	G1531-60630	1
7B. TCD Manifold Regulator Brazement	G1532-80530	1
8. TCD 3-Way Switching Valve	G1532-60570	1




Figure 1320-8 (TCD) Manifold Assembly