

## Capillary Injection Port (0–30 PSIG)

Description	HP Part No.	Qty.
1. SCR-Mach M4 x 0.7 x 8 mm LG	0515-0106	2
2. SCR-Mach M4 x 12 mm, T20 Torx	1390-1023	3
3. Thread Cutting Screw	0624-0665	2

---

**Capillary Injection Port Sub-Assembly (0-30 PSIG)**


---

**Capillary Injection Port Sub-Assembly (0-30 PSIG)**  
**(HP part no. G1544-65500)**

Description	HP Part No.	Qty.
1. O-Ring .239" (Order No. 5180-4182 12/pkg)	0905-1014	1
2. SCR-M4 x 8 T20 Torx-Chrm Plt	0515-2711	2
3. Reducing Nut	18740-20800	1
4. Column Nut	5181-8830	2/pk
5. Gold Plated Seal	18740-20885	1
6. Septum Nut Assembly	18740-60835	1
7. Vent Torque Plate	19244-00030	1
8. Purge Regulator	19246-60530	1
9. Insert Assembly Series II	19251-60575	1
10. Washer, Flat SS	5061-5869	12/pk
11. 11-mm Solsept (Order No. 5182-0739)		50/pkg
12. Ferrule, No Hole	5181-7458	1
13. Cap Top Insulation	G1544-00010	1
14. Flow Controller Inlet Fitting-Replacement	0100-1595	1
15. Cap Inlet Block Insulation	G1544-00020	1
16. Cap Bottom Insulation	G1544-00030	1
17. Cap Heated Block	G1544-20570	1
18. Cap Thread Insert	G1544-20580	1
19. Cap Thermal Nut	G1544-20590	1
20. S/S Solenoid Cable Assembly	G1544-60520	1
21. Cap Solenoid Valve	G1544-60560	1
22. Cap Inj. Heater/Sensor Assembly	G1544-61140	1
23. Cap Shell Weldment	G1544-80570	1
NS Replacement flow controller O-ring kit	0905-1251	-

See also the parts list for the Split/Splitless Inlet with EPC.