

## **Packed Column Injection Port Supplies**

<b>Description</b>	<b>HP Part No.</b>	<b>Qty.</b>
<b>Purged Packed Column Injection Port Kit</b>	5181-8837	
Contains:		
11-mm septa; Graphite/Vespel ferrules; Viton O-rings		
Deactivated disposable glass inserts		

Description	HP Part No.	Qty.
<b>Septum-Purged Packed Column Port</b>		
1. Septum Nut	18740-60835	
Headspace retainer nut	18740-60830	
Merlin Microseal	5181-8816	
Nonpurging septum nut assembly for manual flow control only	19243-60570	
2. Septa, 11-mm OD, no-hole	5181-1263	25/pkg
3. Top Insert Weldment	19243-80570	
4. Viton O-Rings	5080-8898	12/pkg
5. Vespel Ferrules	5080-8774	10/pkg
6. Nut, 1/4" Brass	5180-4105	10/pkg
7. Liners Without Glass Inserts:		
1/8" column liner	19243-80510	
1/4" column liner	19243-80520	
Liners With Glass Inserts:		
HP Series 530- $\mu$ column liner	19244-80540	
1/8" column liner	19243-80530	
1/4" column liner	19243-80540	
Disposable Glass Inserts	5080-8732	25/pkg
Disposable Glass Inserts (deactivated)	5181-3382	5/pkg