

CLP Pesticides

Column: DB-5
15 m x 0.32 mm I.D., 0.25 μ m
10 m DuraGuard Column

Carrier: Helium at 37.2 cm/sec (measured at 50°C)
50°C for 1 min
50-150°C at 25°/min
150-230°C at 4°/min

Injector: Splitless, 250°C
0.5 min purge activation time

Detector: FID, 300°C

- | | |
|-------------------------|------------------------|
| 1. α -BHC | 11. Dieldrin |
| 2. β -BHC | 12. <i>p,p'</i> -DDE |
| 3. γ -BHC | 13. Endrin |
| 4. δ -BHC | 14. Endosulfan II |
| 5. Heptachlor | 15. <i>p,p'</i> -DDT |
| 6. Aldrin | 16. Endrin aldehyde |
| 7. Heptachlor epoxide | 17. Endosulfan sulfate |
| 8. γ -Chlordane | 18. <i>p,p'</i> -DDT |
| 9. Endosulfan I | 19. Endrin ketone |
| 10. α -Chlordane | 20. Methoxychlor |

