

Alltech® Allsphere™ Columns

Try Alltech's lower cost Allsphere™ columns as a replacement for Spherisorb® columns. If you are not satisfied with Allsphere™ as a direct replacement for your Spherisorb® applications, we will refund your money.

Allsphere™ HPLC Columns

Packing	i.d. x Length	Standard Part No.	Waters® Fittings Part No.
ODS-1, 5µm	4.6 x 150mm	778441	778443
	4.6 x 250mm	778364	778357
ODS-2, 3µm	4.6 x 100mm	778487	778489
	4.6 x 150mm	778558	778560
ODS-2, 5µm	4.6 x 150mm	778545	778547
	4.6 x 250mm	778736	778738
C8, 3µm	4.6 x 100mm	778561	778563
	4.6 x 150mm	778609	778611
C8, 5µm	4.6 x 150mm	778571	778573
	4.6 x 250mm	778744	778746
C6, 5µm	4.6 x 150mm	778211	778217
	4.6 x 250mm	778223	778224
Phenyl, 5µm	4.6 x 150mm	778689	778691
	4.6 x 250mm	778752	778754
Cyano, 5µm	4.6 x 150mm	778713	778715
	4.6 x 250mm	778361	778362
Amino, 5µm	4.6 x 150mm	778739	778741
	4.6 x 250mm	778371	778372
Silica, 3µm	4.6 x 100mm	778382	778383
	4.6 x 150mm	778554	778556
Silica, 5µm	4.6 x 150mm	778389	778391
	4.6 x 250mm	778376	778377
SAX, 5µm	4.6 x 250mm	778765	778767
SCX, 5µm	4.6 x 250mm	778188	778189

Allsphere™ All-Guard™ Cartridges*

Packing	i.d. x Length	Qty.	Part No.
ODS-1, 5µm	4.6 x 7.5mm	3	96402
ODS-2, 5µm	4.6 x 7.5mm	3	96403
C8, 5µm	4.6 x 7.5mm	3	96404
C6, 5µm	4.6 x 7.5mm	3	96405
Phenyl, 5µm	4.6 x 7.5mm	3	96406
Cyano, 5µm	4.6 x 7.5mm	3	96408
Amino, 5µm	4.6 x 7.5mm	3	96409
Silica, 5µm	4.6 x 7.5mm	3	96401
SAX, 5µm	4.6 x 7.5mm	3	96410
SCX, 5µm	4.6 x 7.5mm	3	96411
All-Guard™ Cartridge Holder (Includes Direct-Connect Column Coupler)		ea	80101

*All-Guard™ holder required. Other particle sizes available.