



Cross Reference Guide CP-Wax Columns

Varian	Agilent	Restek	Supelco	Alltech	SGE	Quadrex	Machery- Nagel
<u>CP-Wax 52CB</u>	DB-Wax HP-Wax Innowax PEG-20M	Rtx-Wax Stabilwax	Supelcowax-10	AT-Wax	BP-20	007-CW	Optima Wax
<u>CP-Wax 57 CB</u>	Unique column						
<u>CP-Wax 58 (FFAP) CB</u>	DB-FFAP/ HP-FFAP	Stabilwax	SP-1000	AT-1000	SP-1000	007-FFAP	Optima FFAP