

## Put Agilent innovation and quality to work for your non-Agilent GCs

Agilent CrossLab is a growing portfolio of Agilent supplies critical to instrument performance and productivity. Agilent CrossLab supplies feature easy-to-use packaging and include:

- Autosampler syringes
- Inlet septa
- Column ferrules
- Autosampler vials
- Inlet liners
- Column nuts
- and closures
- Liner O-rings
- And more

All Agilent CrossLab supplies are backed by a 90-day guarantee.

So you can put them to the test risk-free.

Plus, we now offer Agilent CrossLab services – a portfolio of service solutions designed to drive productivity and efficiency across your lab:

- Multi-vendor repair
- Compliance
- Laboratory relocation services
- and maintenance

Information, descriptions, and specifications in this publication are subject to change without notice.  
© Agilent Technologies, Inc. 2012  
Printed in the USA March 2, 2012  
5990-9425EN

Agilent Technologies



Learn more: [www.agilent.com/chem/CrossLab](http://www.agilent.com/chem/CrossLab)  
Buy online: [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)  
Find an Agilent Authorized Distributor in your country:  
[www.agilent.com/chem/wheretobuy](http://www.agilent.com/chem/wheretobuy)



Agilent CrossLab works with

BRUKER | PERKINELMER | SHIMADZU | THERMO SCIENTIFIC | AND MORE



# Agilent CrossLab

Supplies for  
**Former Varian**  
GC Systems

**Quick Reference Guide**




The Measure of Confidence



Agilent Technologies

# Inlet Liners

| Liners for 1177 Split/Splitless Injector Ports                                                                  |         |         |             |             |      |                                  |                               |                                  |                                  |  |
|-----------------------------------------------------------------------------------------------------------------|---------|---------|-------------|-------------|------|----------------------------------|-------------------------------|----------------------------------|----------------------------------|--|
| Description                                                                                                     | ID (mm) | OD (mm) | Length (mm) | Volume (µL) | Unit | Similar to Manufacturer Part No. | Agilent Original Deactivation | Similar to Manufacturer Part No. | Agilent Ultra Inert Deactivation |  |
| <b>Split/Splitless Liners</b>                                                                                   |         |         |             |             |      |                                  |                               |                                  |                                  |  |
|  Single taper                  | 4.0     | 6.3     | 78.5        | 1000        | 5/pk | SG092017                         | 8004-0101                     | RT207992145<br>SG092017          | 8004-0151                        |  |
|  Single taper, with wool       | 4.0     | 6.3     | 78.5        | 1000        | 5/pk | SG092019                         | 8004-0102                     | SG092019                         | 8004-0152                        |  |
|  Double taper                  | 4.0     | 6.3     | 78.5        | 1000        | 5/pk | SG092018                         | 8004-0105                     | SG092018                         | 8004-0155                        |  |
|  Gooseneck, with wool          | 4.0     | 6.5     | 78.5        | 1000        | 5/pk | 392611936                        | 8004-0114                     | 392611936                        | 8004-0170                        |  |
|  Recessed gooseneck, with wool | 4.0     | 6.3     | 78.5        | 1000        | 5/pk | SG092010                         | 8004-0103                     | SG092010                         | 8004-0153                        |  |
|  Gooseneck                     | 2.0     | 6.5     | 78.5        | 250         | 5/pk | 392611926                        | 8004-0119                     | 392611926                        | 8004-0178                        |  |
| <b>Splitless Liners</b>                                                                                         |         |         |             |             |      |                                  |                               |                                  |                                  |  |
|  Straight, with wool           | 4.0     | 6.5     | 78.5        | 1000        | 5/pk | 392611937                        | 8004-0116                     | 392611937                        | 8004-0173                        |  |
|  Gooseneck                     | 4.0     | 6.5     | 78.5        | 1000        | 5/pk | 392611927                        | 8004-0113                     | 392611927                        | 8004-0165                        |  |
| <b>Split Liners</b>                                                                                             |         |         |             |             |      |                                  |                               |                                  |                                  |  |
|  Straight-through              | 4.0     | 6.3     | 78.5        | 1000        | 5/pk | SG092007                         | 8004-0156                     | RT207732145<br>SG092007          | 8004-0156                        |  |
|  Straight, with wool           | 4.0     | 6.3     | 78.5        | 1000        | 5/pk | SG092001<br>392611934            | 8004-0104                     | SG092001<br>392611934            | 8004-0154                        |  |
|  With frit, gooseneck          | 4.0     | 6.3     | 78.5        | 1000        | 5/pk |                                  |                               | RT210462145                      | 8004-0158                        |  |
| <b>Direct Liners</b>                                                                                            |         |         |             |             |      |                                  |                               |                                  |                                  |  |
|  Straight-through              | 1.2     | 6.3     | 78.5        | 90          | 5/pk | SG092016                         | 8004-0107                     | SG092016                         | 8004-0157                        |  |

| Liners for 1078/1079 Injector Ports                                                                      |         |         |             |             |      |                                  |                               |                                  |                                  |  |
|----------------------------------------------------------------------------------------------------------|---------|---------|-------------|-------------|------|----------------------------------|-------------------------------|----------------------------------|----------------------------------|--|
| Description                                                                                              | ID (mm) | OD (mm) | Length (mm) | Volume (µL) | Unit | Similar to Manufacturer Part No. | Agilent Original Deactivation | Similar to Manufacturer Part No. | Agilent Ultra Inert Deactivation |  |
| <b>Split/Splitless Liners</b>                                                                            |         |         |             |             |      |                                  |                               |                                  |                                  |  |
|  Single taper         | 3.4     | 5.0     | 54          | 500         | 5/pk | SG092038                         | 8004-0108                     | RT209012145<br>SG092038          | 8004-0160                        |  |
|  Gooseneck, with wool | 2.0     | 5.0     | 54          | 250         | 5/pk | 392611953                        | 8004-0118                     |                                  |                                  |  |
| <b>Splitless Liners</b>                                                                                  |         |         |             |             |      |                                  |                               |                                  |                                  |  |
|  Single taper         | 2.0     | 5.0     | 54          | 170         | 5/pk | SG092039                         | 8004-0109                     | RT207122145<br>SG092039          | 8004-0161                        |  |
| <b>Split Liners</b>                                                                                      |         |         |             |             |      |                                  |                               |                                  |                                  |  |
|  Gooseneck            | 3.4     | 5.0     | 54          | 500         | 5/pk | 392611945                        | 8004-0112                     | 392611945                        | 8004-0164                        |  |
|  With frit, gooseneck | 3.4     | 5.0     | 54          | 500         | 5/pk |                                  |                               | RT217092145                      | 8004-0159                        |  |
|  With frit, gooseneck | 3.4     | 5.0     | 54          | 500         | 5/pk |                                  |                               | 392611946                        | 8004-0171                        |  |
| <b>Other Liners</b>                                                                                      |         |         |             |             |      |                                  |                               |                                  |                                  |  |
|  SPME, straight       | 0.8     | 5.0     | 54          | 30          | 5/pk |                                  |                               | 392611948                        | 8004-0176                        |  |

Varian GC systems are now Bruker products

## Varian GC systems are now Bruker products

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Inlet Liners (Continued)

| Liners for 1093/1094 Injector Ports                                                                                                    |         |         |             |             |      |                                  |                               |                                  |                                  |
|----------------------------------------------------------------------------------------------------------------------------------------|---------|---------|-------------|-------------|------|----------------------------------|-------------------------------|----------------------------------|----------------------------------|
| Description                                                                                                                            | ID (mm) | OD (mm) | Length (mm) | Volume (µL) | Unit | Similar to Manufacturer Part No. | Agilent Original Deactivation | Similar to Manufacturer Part No. | Agilent Ultra Inert Deactivation |
| <b>Direct Liners</b>                                                                                                                   |         |         |             |             |      |                                  |                               |                                  |                                  |
|  SPI for 0.25/0.32 mm id columns                      | 0.5     | 4.6     | 54          | 10          | 5/pk |                                  |                               | 190010906                        | 8004-0167                        |
|  SPI with 0.5 mm restriction for 0.53 mm id on-column | 0.8     | 4.6     | 54          | 30          | 5/pk | SG092034<br>190010907            | 8004-0110                     | SG092034<br>190010907            | 8004-0162                        |

| Liners for 1075/1077 Injector Ports                                                         |         |         |             |             |      |                                  |                               |                                  |                                  |
|---------------------------------------------------------------------------------------------|---------|---------|-------------|-------------|------|----------------------------------|-------------------------------|----------------------------------|----------------------------------|
| Description                                                                                 | ID (mm) | OD (mm) | Length (mm) | Volume (µL) | Unit | Similar to Manufacturer Part No. | Agilent Original Deactivation | Similar to Manufacturer Part No. | Agilent Ultra Inert Deactivation |
| <b>Split Liners</b>                                                                         |         |         |             |             |      |                                  |                               |                                  |                                  |
|  With wool | 4.0     | 6.3     | 72          | 1000        | 5/pk | SG092021<br>190010901            | 8004-0111                     | SG092021<br>190010901            | 8004-0163                        |

| Liners for 1060/1061 Injector Ports                                                                |         |         |             |             |      |                                  |                               |                                  |                                  |
|----------------------------------------------------------------------------------------------------|---------|---------|-------------|-------------|------|----------------------------------|-------------------------------|----------------------------------|----------------------------------|
| Description                                                                                        | ID (mm) | OD (mm) | Length (mm) | Volume (µL) | Unit | Similar to Manufacturer Part No. | Agilent Original Deactivation | Similar to Manufacturer Part No. | Agilent Ultra Inert Deactivation |
| <b>Direct Liners</b>                                                                               |         |         |             |             |      |                                  |                               |                                  |                                  |
|  Double gooseneck | 0.9     | 6.3     | 72          | 1000        | 5/pk |                                  |                               | 392611943                        | 8004-0168                        |

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.



### Ultra Inert deactivation for improved performance and confidence. Every Agilent CrossLab Ultra Inert Liner.

- Achieves accurate analysis of active compounds at nanogram levels, even with glass wool.
- Increases lab productivity by reducing maintenance and recalibration.
- Is certified for consistent quality and reproducible results.

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

# Liner O-Rings

| Liner O-rings                                                         |       |                                  |                           |
|-----------------------------------------------------------------------|-------|----------------------------------|---------------------------|
| Description                                                           | Unit  | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
| Non-stick fluoroelastomer O-ring, 1177 split/splitless, 6.3/6.5 mm od | 10/pk | 8850103100                       | 8004-0201                 |
| Graphite O-ring, 1177 split/splitless, 6.5 mm od                      | 10/pk | 392611930                        | 8004-0202                 |
| Graphite O-ring, 1177 split/splitless, 6.3 mm od                      | 10/pk | 392611935                        | 8004-0203                 |
| Graphite liner seal, 1078/1079 injector, 5 mm id                      | 10/pk | 392534201                        | 8004-0204                 |



Graphite liner O-ring, 8004-0202

# Capillary Column Ferrules

| Injector                                                 | Fitting Size (in.) | Ferrule ID (mm) | Column ID (mm)  | Hole | Unit  | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
|----------------------------------------------------------|--------------------|-----------------|-----------------|------|-------|----------------------------------|---------------------------|
| <b>60% Vespel/40% Graphite Capillary Column Ferrules</b> |                    |                 |                 |      |       |                                  |                           |
| 1177, 1079                                               | 1/16               | 0.3             | 0.18 or smaller | 1    | 10/pk | CR213103                         | 8004-0211                 |
|                                                          | 1/16               | 0.4             | 0.25            | 2    | 10/pk | CR213124                         | 8004-0213                 |
|                                                          | 1/16               | 0.425           | 0.25            | 1    | 10/pk | CR213104                         | 8004-0212                 |
|                                                          | 1/16               | 0.5             | 0.32            | 1    | 10/pk | CR213105                         | 8004-0214                 |
|                                                          | 1/16               | 0.5             | 0.32            | 2    | 10/pk | CR213125                         | 8004-0215                 |
| 1177, 1079, 1061, 1041                                   | 1/16               | 0.8             | 0.53            | 1    | 10/pk | CR213108                         | 8004-0216                 |
| <b>Vespel Capillary Column Ferrules</b>                  |                    |                 |                 |      |       |                                  |                           |
| 1177, 1079                                               | 1/16               | 0.3             | 0.18 or smaller | 1    | 10/pk | CR212103                         | 8010-0306                 |
|                                                          | 1/16               | 0.4             | 0.25            | 1    | 10/pk |                                  | 8010-0307                 |
|                                                          | 1/16               | 0.425           | 0.25            | 1    | 10/pk | CR212104                         | 8004-0219                 |
|                                                          | 1/16               | 0.5             | 0.32            | 1    | 10/pk | CR212105                         | 8010-0308                 |
|                                                          | 1/16               | 0.5             | 0.32            | 2    | 10/pk | CR212125                         | 8004-0218*                |
| 1177, 1079, 1061, 1041                                   | 1/16               | 0.8             | 0.53            | 1    | 10/pk | CR212108                         | 8010-0309                 |
| <b>Graphite Capillary Column Ferrules</b>                |                    |                 |                 |      |       |                                  |                           |
| 1177, 1079                                               | 1/16               | 0.4             | 0.25            | 1    | 10/pk | CR211104                         | 8010-0301                 |
|                                                          | 1/16               | 0.5             | 0.32            | 1    | 10/pk | CR211105                         | 8010-0302                 |
|                                                          | 1/16               | 0.5             | 0.32            | 2    | 10/pk | CR211125                         | 8010-0303                 |
| 1077, 1079, 1061, 1041                                   | 1/16               | 0.8             | 0.53            | 1    | 10/pk | CR211108                         | 8010-0304                 |

\*1177 Injector only

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Packed Column Ferrules

| Injector                                              | Fitting Size (in.) | Ferrule ID (in.) | Column OD (in.) | Hole | Unit  | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
|-------------------------------------------------------|--------------------|------------------|-----------------|------|-------|----------------------------------|---------------------------|
| <b>60% Vespel/40% Graphite Packed Column Ferrules</b> |                    |                  |                 |      |       |                                  |                           |
| 1093, 1061,1041                                       | 1/4                | 1/4              | 1/4             | 1    | 10/pk | CR213400                         | 8004-0217*                |
| <b>Graphite Packed Column Ferrules</b>                |                    |                  |                 |      |       |                                  |                           |
| 1093, 1061,1041                                       | 1/4                | 1/4              | 1/4             | 1    | 10/pk | CR211400                         | 8010-0305*                |

\*Straight body

## Column Nuts

| <b>Column Nuts</b>                                    |      |                                  |                           |
|-------------------------------------------------------|------|----------------------------------|---------------------------|
| Description                                           | Unit | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
| Column nut, brass, 1177, 1079, 1061, or 1041 injector | 2/pk | 394955100                        | 8004-0311                 |
| Column nut, stainless steel, 1093 injector            | 2/pk | CP743117                         | 8004-0312                 |

## Autosampler Syringes

| <b>Autosampler Syringes</b>            |             |                             |                              |      |                                  |                                   |                                              |                                               |
|----------------------------------------|-------------|-----------------------------|------------------------------|------|----------------------------------|-----------------------------------|----------------------------------------------|-----------------------------------------------|
| Model                                  | Volume (µL) | Description                 | Needle Gauge/Length (mm)/Tip | Unit | Similar to Mfr. Syringe Part No. | Agilent CrossLab Syringe Part No. | Agilent CrossLab Replacement Needle Part No. | Agilent CrossLab Replacement Plunger Part No. |
| Varian* CP8400, CP8410, CP9010, CP9050 | 10          | Fixed needle                | 23/50/cone tip               | 1/pk | SG002981                         | 8010-0351                         |                                              |                                               |
|                                        | 10          | Fixed needle                | 26/50/bevel tip              | 1/pk |                                  | 8004-0001                         |                                              |                                               |
|                                        | 10          | Removable needle            | 26/50/cone tip               | 1/pk | SG002982                         | 8004-0003                         | 8004-0004, 2/pk                              |                                               |
| Varian* 8035, 8100, 8200               | 10          | Fixed needle, gas tight     | 26/53/side hole tip          | 1/pk |                                  | 8004-0002                         |                                              | 8004-0007                                     |
|                                        | 10          | Removable needle, gas tight | 25/53/side hole tip          | 1/pk |                                  | 8004-0006                         | 8004-0006, 2/pk                              | 8004-0007                                     |

\*Varian GC systems are now Bruker products

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.



Autosampler syringe, 10 µL fixed needle, 8010-0351

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

# Inlet Septa

| Description                                                  | Similar to Manufacturer Part No. | Agilent CrossLab Part No. 50/pk | Agilent CrossLab Part No. 100/pk |
|--------------------------------------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <b>Non-Stick Bleed and Temperature Optimized (BTO) Septa</b> |                                  |                                 |                                  |
| 9 mm, CenterGuide                                            | CR298713                         | 8010-0217                       | 8010-0218                        |
| 9.5 mm                                                       | CR298705                         | 8010-0219                       | 8010-0220                        |
| 10 mm                                                        | CR298745                         | 8010-0221                       | 8010-0222                        |
| 11 mm, CenterGuide                                           | CR298717                         | 8010-0223                       | 8010-0224                        |
| 11.5 mm, CenterGuide                                         | CR298777                         | 8010-0225                       | 8010-0226                        |
| <b>Non-Stick Advanced Green Septa</b>                        |                                  |                                 |                                  |
| 9 mm, CenterGuide                                            | CR246713                         | 8010-0201                       | 8010-0202                        |
| 9.5 mm                                                       | CR246124                         | 8010-0203                       | 8010-0204                        |
| 10 mm                                                        |                                  | 8010-0205                       | 8010-0206                        |
| 11 mm, CenterGuide                                           | CR246225                         | 8010-0207                       | 8010-0208                        |
| 11.5 mm, CenterGuide                                         | CR246725                         | 8010-0209                       | 8010-0210                        |
| <b>Non-Stick Long-Life Septa</b>                             |                                  |                                 |                                  |
| 9 mm, CenterGuide                                            | CR239778                         | 8010-0233                       | 8010-0234                        |
| 11 mm, CenterGuide                                           | CR239287                         | 8010-0239                       | 8010-0240                        |
| 11.5 mm, CenterGuide                                         | CR239787                         | 8010-0241                       | 8010-0242                        |
| <b>General Purpose Septa</b>                                 |                                  |                                 |                                  |
| 9 mm                                                         |                                  | 8010-0249                       | 8010-0250                        |
| 9.5 mm                                                       |                                  | 8010-0251                       | 8010-0252                        |
| 10 mm                                                        |                                  | 8010-0253                       | 8010-0254                        |
| 11 mm                                                        |                                  | 8010-0255                       | 8010-0256                        |
| 11.5 mm                                                      |                                  | 8010-0257                       | 8010-0258                        |



Non-stick advanced green septum, 11 mm, CenterGuide, 8010-0207

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.



## Innovative, easy-to-use packaging

Agilent CrossLab GC supplies improve overall lab productivity and throughput with:

- Individual blister packs for non-stick plasma-treated inlet septa.
- Touchless packaging of inlet liners.
- Individual dial packs for liner O-rings and column ferrules.
- Proprietary crush barrier for vials.

Pictured above: Non-stick bleed and temperature optimized septa, 10 mm, 50/pk, 8010-0221

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

# Vials and Closures

| 2 mL Wide Opening (9 mm) Screw Top Glass Vials |         |                                       |                           |
|------------------------------------------------|---------|---------------------------------------|---------------------------|
| Description                                    | Unit    | Similar to Manufacturer Part No.      | Agilent CrossLab Part No. |
| Clear                                          | 100/pk  | 392611646<br>392611653                | 8010-0015                 |
| Clear, write-on spot                           | 100/pk  | 39260500<br>CPLC40111B<br>CPLC40111R  | 8010-0014                 |
|                                                | 1000/cs | 392620500<br>CPLC40111B<br>CPLC40111R | 8010-0175                 |
| Amber                                          | 100/pk  | 392611654                             | 8010-0017                 |
| Amber, write-on spot                           | 100/pk  |                                       | 8010-0016                 |
|                                                | 1000/cs |                                       | 8010-0176                 |
| Screw Top Vials with Fixed Inserts             |         |                                       |                           |
| Clear, 300 µL insert volume                    | 100/pk  | 392620514                             | 8010-0008                 |
| Amber, 300 µL insert volume                    | 100/pk  |                                       | 8010-0009                 |

| Screw Caps for Wide Opening (9 mm) Vials |                        |         |                                  |                           |
|------------------------------------------|------------------------|---------|----------------------------------|---------------------------|
| Color                                    | Septa Type             | Unit    | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
| Blue                                     | PTFE/silicone          | 100/pk  | 392611659                        | 8010-0086                 |
|                                          | PTFE/silicone/PTFE     | 100/pk  |                                  | 8010-0087                 |
|                                          | Pre-slit PTFE/silicone | 100/pk  |                                  | 8010-0078                 |
|                                          | Bonded PTFE/silicone   | 100/pk  |                                  | 8010-0077                 |
|                                          | PTFE lined             | 100/pk  |                                  | 8010-0080                 |
|                                          | Open top, no septa     | 100/pk  |                                  | 8010-0079                 |
| Red                                      | PTFE/silicone          | 100/pk  | 392611659                        | 8010-0088                 |
|                                          | PTFE/silicone/PTFE     | 100/pk  |                                  | 8010-0089                 |
|                                          | Pre-slit PTFE/silicone | 100/pk  |                                  | 8010-0082                 |
|                                          | Bonded PTFE/silicone   | 100/pk  |                                  | 8010-0081                 |
|                                          | Open top, no septa     | 100/pk  |                                  | 8010-0083                 |
| Orange                                   | PTFE/silicone          | 100/pk  | 392611659                        | 8010-0084                 |
|                                          |                        | 1000/cs | 392620506                        | 8010-0186                 |
|                                          | PTFE/silicone/PTFE     | 100/pk  |                                  | 8010-0085                 |
|                                          |                        | 1000/cs |                                  | 8010-0187                 |
|                                          | Pre-slit PTFE/silicone | 100/pk  |                                  | 8010-0075                 |
|                                          | Bonded PTFE/silicone   | 100/pk  |                                  | 8010-0074                 |
|                                          | Open top, no septa     | 100/pk  |                                  | 8010-0076                 |

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

## NOTE:

To obtain an optimal fit and seal of any vial and closure, we strongly recommend buying an entire Agilent CrossLab vial assembly (including a vial, cap, and septum, as well as a vial insert if needed).

For additional help in selecting vials, caps, and septa, go to [www.agilent.com/chem/SelectVials](http://www.agilent.com/chem/SelectVials)



Clear, 2 mL wide opening (9 mm) screw top glass vials, 8010-0015



Orange screw caps with PTFE/silicone septa for wide opening (9 mm) vials, 8010-0084

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Vials and Closures *(Continued)*

### Septa for Wide Opening (9 mm) Screw Caps

| Septa Type                   | Color | Unit    | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
|------------------------------|-------|---------|----------------------------------|---------------------------|
| PTFE/red silicone            | Ivory | 100/pk  | 190010209<br>392620510           | 8010-0093                 |
| PTFE/white silicone/red PTFE | Red   | 100/pk  |                                  | 8010-0091                 |
|                              |       | 1000/cs |                                  | 8010-0188                 |
| Pre-slit PTFE/white silicone | Blue  | 100/pk  |                                  | 8010-0094                 |
| PTFE/white silicone          | Red   | 100/pk  | 190010209<br>392620510           | 8010-0092                 |

### 2 mL Wide Opening (9 mm) Screw Top Vial Convenience Packs

| Vial Type            | Septa Type    | Cap Color | Unit   | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
|----------------------|---------------|-----------|--------|----------------------------------|---------------------------|
| Clear, write-on spot | PTFE/silicone | Orange    | 100/pk | 392611549                        | 8010-0198                 |
| Amber, write-on spot | PTFE/silicone | Orange    | 100/pk | 392611550                        | 8010-0199                 |

### 2 mL Standard Opening (8 mm) Screw Top Glass Vials

| Description          | Unit   | Similar to Manufacturer Part No.  | Agilent CrossLab Part No. |
|----------------------|--------|-----------------------------------|---------------------------|
| Clear                | 100/pk | 392611639<br>392620550<br>CP10271 | 8010-0011                 |
| Clear, write-on spot | 100/pk |                                   | 8010-0010                 |
| Amber                | 100/pk | 392611640                         | 8010-0013                 |
| Amber, write-on spot | 100/pk |                                   | 8010-0012                 |

### 2 mL Standard Opening (8 mm) Screw Top Glass Vial Convenience Packs

| Vial Type | Septa Type    | Cap Color | Unit   | Similar to Manufacturer Part No.                 | Agilent CrossLab Part No. |
|-----------|---------------|-----------|--------|--------------------------------------------------|---------------------------|
| Clear     | PTFE/silicone | Black     | 100/pk | 392611530<br>392611987<br>394983500<br>394983501 | 8010-0414                 |
| Amber     | PTFE/silicone | Black     | 100/pk | 392611531                                        | 8010-0415                 |

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are a viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

#### NOTE:

To obtain an optimal fit and seal of any vial and closure, we strongly recommend buying an entire Agilent CrossLab vial assembly (including a vial, cap, and septum, as well as a vial insert if needed).

For additional help in selecting vials, caps, and septa, go to [www.agilent.com/chem/SelectVials](http://www.agilent.com/chem/SelectVials)



Standard opening screw top vials, 8010-0011

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Vials and Closures (Continued)

| Screw Caps for Standard Opening (8 mm) Vials |                                 |        |                                  |                           |
|----------------------------------------------|---------------------------------|--------|----------------------------------|---------------------------|
| Color                                        | Septa Type                      | Unit   | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
| Blue                                         | PTFE/silicone                   | 100/pk | 392611641                        | 8010-0061                 |
|                                              | Pre-slit PTFE/silicone          | 100/pk |                                  | 8010-0062                 |
|                                              | PTFE/silicone/PTFE              | 100/pk |                                  | 8010-0063                 |
|                                              | PTFE/butyl                      | 100/pk | 392611645                        | 8010-0064                 |
|                                              | Open top, no septa              | 100/pk |                                  | 8010-0065                 |
| Red                                          | PTFE/silicone                   | 100/pk | 392611641                        | 8010-0066                 |
|                                              | Pre-slit PTFE/silicone          | 100/pk |                                  | 8010-0067                 |
|                                              | PTFE/silicone/PTFE              | 100/pk |                                  | 8010-0068                 |
|                                              | PTFE/butyl                      | 100/pk | 392611645                        | 8010-0069                 |
|                                              | Open top, no septa              | 100/pk | 1600069800                       | 8010-0070                 |
| Orange                                       | PTFE/silicone                   | 100/pk | 392611641                        | 8010-0056                 |
|                                              | Pre-slit PTFE/silicone          | 100/pk |                                  | 8010-0057                 |
|                                              | PTFE/silicone/PTFE              | 100/pk |                                  | 8010-0058                 |
|                                              | PTFE/butyl                      | 100/pk | 392611645                        | 8010-0059                 |
|                                              | Open top, no septa              | 100/pk |                                  | 8010-0060                 |
| Black                                        | PTFE/silicone                   | 100/pk |                                  | 8010-0054                 |
|                                              | Open top, no septa              | 100/pk |                                  | 8010-0055                 |
| Flanged Caps                                 |                                 |        |                                  |                           |
| Orange                                       | Flanged, PTFE/silicone          | 100/pk |                                  | 8010-0141                 |
|                                              | Flanged, pre-slit PTFE/silicone | 100/pk | 392611613                        | 8010-0142                 |
|                                              | Flanged, no septa               | 100/pk |                                  | 8010-0166                 |

| Septa for Standard Opening (8 mm) Vials |        |                                       |                           |
|-----------------------------------------|--------|---------------------------------------|---------------------------|
| Septa Type                              | Unit   | Similar to Manufacturer Part No.      | Agilent CrossLab Part No. |
| Red PTFE/white silicone/red PTFE        | 100/pk |                                       | 8010-0071                 |
| Red PTFE/white silicone                 | 100/pk | 390614901<br>6900016900<br>R005486RT1 | 8010-0072                 |
| PTFE/butyl                              | 500/pk |                                       | 8010-0073                 |

| 4 mL (13 mm) Screw Top Glass Vials |        |                                  |                           |
|------------------------------------|--------|----------------------------------|---------------------------|
| Description                        | Unit   | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
| Clear                              | 100/pk | CP10277                          | 8010-0020                 |
| Amber                              | 100/pk |                                  | 8010-0021                 |

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

### NOTE:

To obtain an optimal fit and seal of any vial and closure, we strongly recommend buying an entire Agilent CrossLab vial assembly (including a vial, cap, and septum, as well as a vial insert if needed).

For additional help in selecting vials, caps, and septa, go to [www.agilent.com/chem/SelectVials](http://www.agilent.com/chem/SelectVials)



Screw caps for standard opening vials, 8010-0056

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Vials and Closures (Continued)

| Screw Caps for 4 mL (13 mm) Vials |                    |        |                           |
|-----------------------------------|--------------------|--------|---------------------------|
| Color                             | Septa Type         | Unit   | Agilent CrossLab Part No. |
| Black                             | Red PTFE/silicone  | 100/pk | 8010-0095                 |
|                                   | Open top, no septa | 100/pk | 8010-0096                 |

| Septa for 4 mL (13 mm) Vials |        |                           |
|------------------------------|--------|---------------------------|
| Septa Type                   | Unit   | Agilent CrossLab Part No. |
| Red PTFE/white silicone      | 100/pk | 8010-0098                 |
| White PTFE                   | 100/pk | 8010-0099                 |

| 6 mL (16 mm) Screw Top Glass Vials |        |                           |
|------------------------------------|--------|---------------------------|
| Description                        | Unit   | Agilent CrossLab Part No. |
| Clear                              | 100/pk | 8010-0022                 |

| Screw Caps for 6 mL (16 mm) Vials |        |                           |
|-----------------------------------|--------|---------------------------|
| Septa Type                        | Unit   | Agilent CrossLab Part No. |
| PTFE/silicone                     | 100/pk | 8010-0101                 |
| Pre-slit PTFE/silicone            | 100/pk | 8010-0102                 |
| Open top, no septa                | 100/pk | 8010-0100                 |

| 2 mL (11 mm) Crimp Top Glass Vials |         |                                  |                           |
|------------------------------------|---------|----------------------------------|---------------------------|
| Description                        | Unit    | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
| Clear                              | 100/pk  | 390886500*<br>392611634          | 8010-0002                 |
| Clear, write-on spot               | 100/pk  | 392620517                        | 8010-0001                 |
|                                    | 1000/cs | 392620517                        | 8010-0170                 |
| Amber, write-on spot               | 1000/cs | CPLCL07CPVA                      | 8010-0172                 |
|                                    | 100/pk  |                                  | 8010-0003                 |
| Crimp Top Vials with Fixed Inserts |         |                                  |                           |
| Clear, 300 µL insert volume        | 100/pk  |                                  | 8010-0006                 |
| Amber, 300 µL insert volume        | 100/pk  |                                  | 8010-0007                 |

\*144/pk

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

### NOTE:

To obtain an optimal fit and seal of any vial and closure, we strongly recommend buying an entire Agilent CrossLab vial assembly (including a vial, cap, and septum, as well as a vial insert if needed).

For additional help in selecting vials, caps, and septa, go to [www.agilent.com/chem/SelectVials](http://www.agilent.com/chem/SelectVials)



Black screw caps with red PTFE/silicone septa for 4 mL (13 mm) vials, 8010-0095



Crimp top glass vials with write-on spot, 8010-0001

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Vials and Closures (Continued)

### Crimp Caps with 11 mm Septa

| Cap Color       | Septa Type           | Unit    | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
|-----------------|----------------------|---------|----------------------------------|---------------------------|
| Silver aluminum | PTFE/silicone/PTFE   | 100/pk  |                                  | 8010-0049                 |
|                 |                      | 1000/cs |                                  | 8010-0182                 |
|                 | PTFE/silicone        | 100/pk  | 392611631                        | 8010-0050                 |
|                 |                      | 1000/cs | CPLCL11ACBST144                  | 8010-0183                 |
|                 | PTFE/butyl           | 100/pk  | 392611632                        | 8010-0051                 |
|                 | Thin membrane rubber | 25/pk   |                                  | 8010-0053                 |
|                 | PTFE/red rubber      | 100/pk  |                                  | 8010-0045                 |
|                 |                      | 1000/cs |                                  | 8010-0181                 |
| Blue aluminum   | PTFE/red rubber      | 100/pk  |                                  | 8010-0046                 |
| Green aluminum  | PTFE/red rubber      | 100/pk  |                                  | 8010-0047                 |
| Red aluminum    | PTFE/red rubber      | 100/pk  |                                  | 8010-0048                 |
| Gold aluminum   | Magnetic             | 100/pk  |                                  | 8010-0052                 |

### Crimp Top Micro Vials

| Description           | Volume | Unit    | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
|-----------------------|--------|---------|----------------------------------|---------------------------|
| Amber, flat bottom    | 0.8 mL | 1000/pk |                                  | 8010-0153                 |
| Clear, tapered bottom | 0.1 mL | 500/pk  |                                  | 8010-0154                 |
| Clear, rounded bottom | 0.3 mL | 500/pk  | CP959485                         | 8010-0155                 |

### 2 mL (11 mm) Crimp Top Vial Convenience Packs

| Vial Type            | Septa Type    | Cap Color       | Unit   | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
|----------------------|---------------|-----------------|--------|----------------------------------|---------------------------|
| Clear, write-on spot | PTFE/silicone | Silver aluminum | 100/pk | 392611525                        | 8010-0195                 |
| Amber, write-on spot | PTFE/silicone | Silver aluminum | 100/pk |                                  | 8010-0196                 |

### 2 mL (11 mm) Snap Top Glass Vials

| Description          | Unit    | Agilent CrossLab Part No. |
|----------------------|---------|---------------------------|
| Clear                | 100/pk  | 8010-0024                 |
| Clear, write-on spot | 100/pk  | 8010-0023                 |
|                      | 1000/cs | 8010-0177                 |
| Amber                | 100/pk  | 8010-0026                 |
| Amber, write-on spot | 100/pk  | 8010-0025                 |
|                      | 1000/cs | 8010-0178                 |

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

#### NOTE:

To obtain an optimal fit and seal of any vial and closure, we strongly recommend buying an entire Agilent CrossLab vial assembly (including a vial, cap, and septum, as well as a vial insert if needed).

For additional help in selecting vials, caps, and septa, go to [www.agilent.com/chem/SelectVials](http://www.agilent.com/chem/SelectVials)



Amber, flat bottom crimp top micro vials, 8010-0153



Amber, write-on spot 2 mL (11 mm) snap top glass vials, 8010-0025

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Vials and Closures (Continued)

| Snap Caps with 11 mm Septa |                              |         |                                  |                           |
|----------------------------|------------------------------|---------|----------------------------------|---------------------------|
| Cap Color                  | Septa Type                   | Unit    | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
| Clear                      | PTFE/red silicone            | 100/pk  |                                  | 8010-0104                 |
|                            | PTFE/silicone/PTFE           | 100/pk  |                                  | 8010-0105                 |
|                            | Pre-slit PTFE/white silicone | 100/pk  |                                  | 8010-0106                 |
|                            | Clear polyethylene membrane  | 100/pk  |                                  | 8010-0107                 |
| Blue                       | PTFE/red silicone            | 100/pk  | 392611871                        | 8010-0108                 |
|                            | PTFE/white silicone          | 100/pk  |                                  | 8010-0112                 |
| Green                      | PTFE/red silicone            | 100/pk  |                                  | 8010-0109                 |
| Orange                     | PTFE/red silicone            | 100/pk  |                                  | 8010-0103                 |
|                            | PTFE/red silicone            | 1000/cs |                                  | 8010-0189                 |
|                            | PTFE/white silicone          | 100/pk  |                                  | 8010-0111                 |
|                            | Pre-slit PTFE/white silicone | 100/pk  |                                  | 8010-0113                 |
|                            | Pre-slit PTFE/white silicone | 1000/cs |                                  | 8010-0190                 |
| Red                        | PTFE/red silicone            | 100/pk  |                                  | 8010-0110                 |
| Polyurethane               |                              | 100/pk  |                                  | 8010-0115                 |

| 2 mL (11 mm) Snap Top Vial Convenience Packs |               |           |        |                           |
|----------------------------------------------|---------------|-----------|--------|---------------------------|
| Vial Type                                    | Septa Type    | Cap Color | Unit   | Agilent CrossLab Part No. |
| Clear, write-on spot                         | PTFE/silicone | Orange    | 100/pk | 8010-0200                 |
| Amber, write-on spot                         | PTFE/silicone | Orange    | 100/pk | 8010-0411                 |

| Snap Top Solvent Vial Convenience Packs      |                         |           |        |                                                                            |                           |
|----------------------------------------------|-------------------------|-----------|--------|----------------------------------------------------------------------------|---------------------------|
| Vial Type                                    | Septa Type              | Cap Color | Unit   | Similar to Manufacturer Part No.                                           | Agilent CrossLab Part No. |
| 5 mL Snap Top Solvent Vial Convenience Pack  |                         |           |        |                                                                            |                           |
| Clear                                        | PTFE/silicone           | Clear     | 100/pk | 392611970<br>392611974<br>392611975<br>392611980<br>392611981<br>392611982 | 8010-0416                 |
| 10 mL Snap Top Solvent Vial Convenience Pack |                         |           |        |                                                                            |                           |
| Clear                                        | Red PTFE/white silicone | Clear     | 100/pk | 392612006<br>392612007<br>392611972<br>392611973                           | 8010-0422                 |

### NOTE:

To obtain an optimal fit and seal of any vial and closure, we strongly recommend buying an entire Agilent CrossLab vial assembly (including a vial, cap, and septum, as well as a vial insert if needed).

For additional help in selecting vials, caps, and septa, go to [www.agilent.com/chem/SelectVials](http://www.agilent.com/chem/SelectVials)



Orange snap caps with PTFE/red silicone 11 mm septa, 8010-0103

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Vials and Closures (Continued)

### Snap Cap for 5 mL and 10 mL Solvent Vials

| Cap Color | Septa Type                     | Unit  | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
|-----------|--------------------------------|-------|----------------------------------|---------------------------|
| Clear     | Red ultra PTFE/yellow silicone | 50/pk | 392611984                        | 8010-0423                 |

### 19 mm Septa for 5 mL and 10 mL Snap Caps

| Septa Type | Unit   | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
|------------|--------|----------------------------------|---------------------------|
| Red rubber | 100/pk | CP10225<br>392611986             | 8010-0430                 |

### Polypropylene Vials

| Description                                     | Recommended Fill Volume | Max Fill Volume | Unit    | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
|-------------------------------------------------|-------------------------|-----------------|---------|----------------------------------|---------------------------|
| Polypropylene, crimp/snap top                   | 250 µL                  | 440 µL          | 100/pk  |                                  | 8010-0158                 |
| Polypropylene, crimp/snap top                   | 250 µL                  | 440 µL          | 1000/cs |                                  | 8010-0193                 |
| Polypropylene, crimp/snap top                   | 0.7 mL                  | 1 mL            | 100/pk  | 392611667                        | 8010-0159                 |
| Polypropylene, screw top                        | 250 µL                  | 440 µL          | 1000/cs |                                  | 8010-0194                 |
| Polypropylene, screw top                        | 250 µL                  | 440 µL          | 100/pk  |                                  | 8010-0161                 |
| Polypropylene with glass insert, crimp/snap top | 250 µL                  | 350 µL          | 100/pk  | R000548SCY                       | 8010-0160                 |
| Polypropylene with glass insert, screw top      | 250 µL                  | 350 µL          | 100/pk  |                                  | 8010-0162                 |

### Micro-V Vials

| Description      | Unit   | Agilent CrossLab Part No. |
|------------------|--------|---------------------------|
| Clear, screw top | 100/pk | 8010-0018                 |
| Clear, crimp top | 100/pk | 8010-0004                 |
| Clear, snap top  | 100/pk | 8010-0027                 |
| Amber, screw top | 100/pk | 8010-0019                 |
| Amber, crimp top | 100/pk | 8010-0005                 |
| Amber, snap top  | 100/pk | 8010-0028                 |

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

### NOTE:

To obtain an optimal fit and seal of any vial and closure, we strongly recommend buying an entire Agilent CrossLab vial assembly (including a vial, cap, and septum, as well as a vial insert if needed).

For additional help in selecting vials, caps, and septa, go to [www.agilent.com/chem/SelectVials](http://www.agilent.com/chem/SelectVials)



Red rubber septa (19 mm) for 5 mL and 10 mL snap caps, 8010-0430



Polypropylene crimp/snap top vial with glass insert, 8010-0160



Clear, crimp top Micro-V vial, 8010-0004

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Vials and Closures (Continued)

| Microvolume Inserts for 8 mm Screw Top Vials |                         |                 |        |                                  |                           |
|----------------------------------------------|-------------------------|-----------------|--------|----------------------------------|---------------------------|
| Description                                  | Recommended Fill Volume | Max Fill Volume | Unit   | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
| Glass conical insert with polymer feet       | 150 µL                  | 200 µL          | 100/pk |                                  | 8010-0124                 |
| Glass pulled-point insert                    | 150 µL                  | 200 µL          | 100/pk | 392611593                        | 8010-0125                 |
| Flat bottom insert                           | 200 µL                  | 220 µL          | 100/pk |                                  | 8010-0126                 |

| Microvolume Inserts for 11 mm Crimp Top or 9 mm Screw Top Vials |                         |                 |        |                                     |                           |
|-----------------------------------------------------------------|-------------------------|-----------------|--------|-------------------------------------|---------------------------|
| Description                                                     | Recommended Fill Volume | Max Fill Volume | Unit   | Similar to Manufacturer Part No.    | Agilent CrossLab Part No. |
| Conical insert with polymer feet                                | 250 µL                  | 330 µL          | 100/pk | 190010204<br>392611594<br>392620538 | 8010-0131                 |
| Pulled-point insert                                             | 250 µL                  | 350 µL          | 100/pk | 392611596<br>392611636              | 8010-0132                 |
| Flat bottom insert                                              | 400 µL                  | 440 µL          | 500/pk |                                     | 8010-0136                 |

| Headspace Crimp Top Glass Vials       |         |                                     |                           |
|---------------------------------------|---------|-------------------------------------|---------------------------|
| Description                           | Unit    | Similar to Manufacturer Part No.    | Agilent CrossLab Part No. |
| 6 mL Headspace Crimp Top Glass Vials  |         |                                     |                           |
| Clear                                 | 100/pk  | 392611857                           | 8010-0151*                |
| 10 mL Headspace Crimp Top Glass Vials |         |                                     |                           |
| Clear, flat bottom                    | 100/pk  | CP10204                             | 8010-0029                 |
| Clear, flat bottom                    | 1000/cs | CP10204                             | 8010-0179                 |
| Clear, rounded bottom                 | 100/pk  | MLA21000                            | 8010-0030                 |
| Amber, flat bottom                    | 100/pk  |                                     | 8010-0031                 |
| Amber, rounded bottom                 | 100/pk  | 392620526                           | 8010-0032                 |
| 20 mL Headspace Crimp Top Glass Vials |         |                                     |                           |
| Clear, flat bottom                    | 100/pk  | CP10205                             | 8010-0033                 |
| Clear, flat bottom                    | 1000/cs | CP10205                             | 8010-0180                 |
| Clear, rounded bottom                 | 100/pk  | 392611675<br>392612021<br>MLA202100 | 8010-0034                 |
| Amber, flat bottom                    | 100/pk  |                                     | 8010-0035                 |
| Amber, rounded bottom                 | 100/pk  | 392620528                           | 8010-0036                 |
| 22 mL Headspace Crimp Top Vials       |         |                                     |                           |
| Clear                                 | 100/pk  | MLA110040M                          | 8010-0152*                |

\*Use 20 mm crimp caps with 6 mL and 22 mL headspace vials

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

### NOTE:

To obtain an optimal fit and seal of any vial and closure, we strongly recommend buying an entire Agilent CrossLab vial assembly (including a vial, cap, and septum, as well as a vial insert if needed).

For additional help in selecting vials, caps, and septa, go to [www.agilent.com/chem/SelectVials](http://www.agilent.com/chem/SelectVials)



Amber, rounded bottom 10 mL headspace crimp top glass vials, 8010-0032

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Vials and Closures (Continued)

| Headspace Screw Top Glass Vials              |        |                                  |                           |
|----------------------------------------------|--------|----------------------------------|---------------------------|
| Description                                  | Unit   | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
| <b>10 mL Headspace Screw Top Glass Vials</b> |        |                                  |                           |
| Clear, flat bottom                           | 100/pk |                                  | 8010-0037                 |
| Clear, rounded bottom                        | 100/pk | CP741315<br>CP913323             | 8010-0038                 |
| Amber, flat bottom                           | 100/pk |                                  | 8010-0039                 |
| Amber, rounded bottom                        | 100/pk | 392620104                        | 8010-0040                 |
| <b>20 mL Headspace Screw Top Glass Vials</b> |        |                                  |                           |
| Clear, flat bottom                           | 100/pk |                                  | 8010-0041                 |
| Clear, rounded bottom                        | 100/pk | 392620202                        | 8010-0042                 |
| Amber, flat bottom                           | 100/pk |                                  | 8010-0043                 |
| Amber, rounded bottom                        | 100/pk | 392620204                        | 8010-0044                 |

| Headspace Caps and Septa                    |                    |         |                                        |                           |
|---------------------------------------------|--------------------|---------|----------------------------------------|---------------------------|
| Cap Color                                   | Septa Type         | Unit    | Similar to Manufacturer Part No.       | Agilent CrossLab Part No. |
| <b>18 mm Headspace Crimp Caps and Septa</b> |                    |         |                                        |                           |
| Silver                                      | Tan PTFE/silicone  | 100/pk  |                                        | 8010-0119                 |
| <b>20 mm Headspace Crimp Caps and Septa</b> |                    |         |                                        |                           |
| Silver aluminum                             | Tan PTFE/silicone  | 100/pk  | 392610003<br>392611858<br>ROK7382420   | 8010-0116                 |
|                                             |                    | 1000/cs | 392610003<br>392611858<br>ROK7382420   | 8010-0191                 |
|                                             | Gray PTFE/butyl    | 100/pk  |                                        | 8010-0117                 |
|                                             | No septa           | 100/pk  |                                        | 8010-0120                 |
| Silver aluminum with safety feature         | PTFE/butyl         | 100/pk  | 392611607                              | 8010-0143                 |
|                                             | PTFE/silicone      | 100/pk  |                                        | 8010-0144                 |
| Silver, magnetic                            | Tan PTFE/silicone  | 100/pk  | 392611860<br>392612017<br>MLA2000051ML | 8010-0165                 |
| Bimetal, magnetic                           | PTFE/silicone      | 100/pk  | CP959464                               | 8010-0420                 |
| Steel                                       | High performance   | 100/pk  |                                        | 8010-0429*                |
| Septa only                                  | Gray PTFE/butyl    | 100/pk  | 392611859                              | 8010-0121                 |
|                                             | Tan PTFE/silicone  | 100/pk  |                                        | 8010-0122                 |
|                                             | Tan PTFE/silicone  | 1000/cs |                                        | 8010-0192                 |
| <b>18 mm Headspace Screw Caps and Septa</b> |                    |         |                                        |                           |
| Silver, magnetic                            | PTFE/silicone      | 100/pk  | 392620302<br>392620304                 | 8010-0139**               |
|                                             | PTFE/butyl         | 100/pk  | 392620306                              | 8010-0140                 |
| Steel                                       | High performance   | 100/pk  |                                        | 8010-0428*                |
| Septa only                                  | Blue PTFE/silicone | 1000/cs | CP912446<br>CP913077                   | 8010-0418                 |

\*Recommended for high temperature applications up to 300 °C for up to one hour.

\*\*This cap is suitable for use in SPME applications.

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

### NOTE:

To obtain an optimal fit and seal of any vial and closure, we strongly recommend buying an entire Agilent CrossLab vial assembly (including a vial, cap, and septum, as well as a vial insert if needed).

For additional help in selecting vials, caps, and septa, go to [www.agilent.com/chem/SelectVials](http://www.agilent.com/chem/SelectVials)



Clear, rounded bottom 20 mL headspace screw top glass vials, 8010-0042



Silver, magnetic 18 mm headspace screw caps with PTFE/butyl septa, 8010-0140

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)

## Vials and Closures (Continued)

| Headspace Vial Convenience Packs        |               |                 |         |                                  |                           |
|-----------------------------------------|---------------|-----------------|---------|----------------------------------|---------------------------|
| Vial Type                               | Septa Type    | Cap Color       | Unit    | Similar to Manufacturer Part No. | Agilent CrossLab Part No. |
| 10 mL, crimp top, clear, flat bottom    | PTFE/silicone | Silver aluminum | 100/pk  | 392611613                        | 8010-0412                 |
| 20 mL, crimp top, clear, flat bottom    | PTFE/silicone | Silver aluminum | 100/pk  | 392611607                        | 8010-0413                 |
| 20 mL, screw top, clear, rounded bottom | PTFE/silicone | Silver magnetic | 1000/cs | 392620401<br>392620403           | 8010-0417                 |

The cross reference to manufacturer part numbers listed here serve as a recommendation that the Agilent CrossLab products are viable alternative to the manufacturer's products. Although Agilent CrossLab products may have, in some cases, slightly different designs as compared to the manufacturer's counterparts, Agilent CrossLab products are compatible with the corresponding Varian (now Bruker) instruments. All Agilent CrossLab supplies are backed by Agilent's 90-day warranty and money-back guarantee.

### NOTE:

To obtain an optimal fit and seal of any vial and closure, we strongly recommend buying an entire Agilent CrossLab vial assembly (including a vial, cap, and septum, as well as a vial insert if needed).

For additional help in selecting vials, caps, and septa, go to [www.agilent.com/chem/SelectVials](http://www.agilent.com/chem/SelectVials)



To learn more about our comprehensive portfolio of supplies for your other non-Agilent GCs, go to [www.agilent.com/chem/CrossLab](http://www.agilent.com/chem/CrossLab)



Or scan now to download our interactive Selection Guide to your mobile device.

### Quality you can count on

Agilent CrossLab supplies are designed to work seamlessly with the non-Agilent instruments in your lab. All Agilent CrossLab supplies are backed by our 90-day warranty and money-back guarantee. That means you have 90 risk-free days to put them to the test!

Contact your Agilent Authorized Distributor, or order now at [www.agilent.com/chem/CrossLabstore](http://www.agilent.com/chem/CrossLabstore)