




YMC-Triart C18



YMC-Triart C18 ExRS




YMC-Triart C8




Easy Column Selection

for Method Development


YMC-Triart Radar Chart



YMC-Triart Phenyl



YMC-Triart PFP



Transfer

Scalable particles:
EASY
UHPLC ↔ HPLC

Flexible

pH = 1 – 12
Temp. up to 70°C
100% aqueous conditions

Universal

YMC-Triart
for acidic, basic and neutral analytes

YMC
EUROPE GMBH

Schöttmannshof 19 · D-46539 Dinslaken
Phone: +49 (0)2064 427-0 · Fax: +49 (0)2064 427-222 · E-Mail: info@ymc.de

www.ymc.de