

# S.C.A.T. Europe Exhaust Filter Overview

## Generation V3.0, Active Carbon 3.0



THE System for Solvent Safety

### Exhaust Filter S

- With change indicator
- Active carbon 3.0
- Economy package available

**3 months\* service life**



• One economy package provides for a year's requirement.

### Exhaust Filter M

- With change indicator
- Active carbon 3.0
- Economy package available

**6 months\* service life**



• One economy package provides for a year's requirement.

### Exhaust Filter L

- With change indicator
- Active carbon 3.0
- Economy package available

**12 months\* service life**



• One economy package provides for 2 years' requirement.

\*With typical HPLC flowrates of 1.5 - 4.0 ml per minute.

Part No.	Filter size	Units	Description
410 534	S	1	Exhaust filter <b>S</b> , V3.0, with splash guard and <b>change indicator</b> , service life <b>3 months</b>
490 335	S	4	Exhaust filter <b>S</b> , V3.0, <b>economy package</b> , with splash guard and <b>change indicator</b> , service life <b>3 months</b>
410 535	M	1	Exhaust filter <b>M</b> , V3.0, with splash guard and <b>change indicator</b> , service life <b>6 months</b>
490 336	M	2	Exhaust filter <b>M</b> , V3.0, <b>economy package</b> , with splash guard and <b>change indicator</b> , service life <b>6 months</b>
407 986	L	1	Exhaust filter <b>L</b> , V3.0, with splash guard and <b>change indicator</b> , service life <b>12 months</b>
490 986	L	2	Exhaust filter <b>L</b> , V3.0, <b>economy package</b> , with splash guard and <b>change indicator</b> , service life <b>12 months</b>

#200web / 02122020