



**THE No. 1 System
for Solvent Safety**

Overview

Range of Canisters & Laboratory Glass Bottles

All-in with SCAT!



SCAT Europe Overview

Range of Canisters



**THE No. 1 System
for Solvent Safety**



Product Group 1: Electrically Conducting Canisters

Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D	Part Number
A	S 50	5 litres	PE-HD-el	black	No	designed as space-saving canister	65 x 335 x 335	108 317
B	S 50	10 litres	PE-HD-el	black	Yes	UN-X approval; with viewing strip	190 x 315 x 230	108 421
C	S 60 / 61	10 litres	PE-HD-el	black	Yes	UN-Y approval	185 x 265 x 290	107 953
D	S 60 / 61	10 litres	PE-HD-el	black	No	with red marking strip	185 x 265 x 290	108 216
E	S 60 / 61	10 litres	PE-HD-el	black	No	with green marking strip	185 x 265 x 290	108 217
F	S 60 / 61	10 litres	PE-HD-el	black	No	with yellow marking strip	185 x 265 x 290	108 215
G	S 60 / 61	10 litres	PE-HD-el	black	No	with blue marking strip	185 x 265 x 290	108 214
H	S 60 / 61	10 litres	PE-HD-el	black	No	with float	185 x 265 x 290	108 042
I	S 60 / 61	10 litres	PE-HD-el	black	No	with holder for sensor 108 178	185 x 265 x 290	108 042-51
J	S 60 / 61	20 litres	PE-HD-el	black	Yes	UN-Y approval	185 x 500 x 290	108 027
K	S 60 / 61	20 litres	PE-HD-el	black	No	with float	185 x 500 x 290	108 043
L	S 60 / 61	30 litres	PE-HD-el	black	Yes	UN-Y approval	240 x 455 x 364	108 192
M	S 60 / 61	30 litres	PE-HD-el	black	No	with float	240 x 455 x 364	109 193
N	S 70 / 71	60 litres	PE-HD-el	black	Yes	UN-Y approval	330 x 625 x 396	107 940
O	S 70 / 71	60 litres	PE-HD-el	black	No	with float	330 x 690 x 395	107 740
P	S 90	10 litres	PE-HD-el	black	Yes	UN-X approval; with viewing strip	195 x 380 x 195	108 420

SCAT Europe Overview

Range of Canisters



**THE No. 1 System
for Solvent Safety**



Product Group 2: Non-Electrically Conducting Canisters



Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D	Part Number
A	GL 45	2,5 litres	PE-HD	natural (white)	Yes	UN-Y approval	115 x 210 x 150	107 950
B	GL 45	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	150 x 250 x 195	107 951
C	GL 45	10 litres	PE-HD	natural (white)	Yes	UN-Y approval	190 x 305 x 230	107 952
D	S 51	5 litres	PE-HD	natural (white)	Yes	UN-X approval	145 x 250 x 190	107 958
E	S 55	5 litres	PE-HD	natural (white)	Yes	UN-X approval	182 x 240 x 162	108 175
F	S 55	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	160 x 230 x 185	107 957
G	S 55	5 litres	F-PE-HD	natural (white)	Yes	UN-X approval; with inner fluorination	145 x 250 x 190	107 933
H	S 55	10 litres	F-PE-HD	natural (white)	Yes	UN-X approval; with inner fluorination	190 x 310 x 230	108 189
I	S 55	10 litres	PE-HD	natural (white)	Yes	UN-Y approval	185 x 305 x 225	107 955
J	S 60 / 61	5 litres	PE-HD	natural (white)	Yes	UN-X approval	164 x 235 x 185	108 227
K	S 60 / 61	12 litres	PE-HD	natural (white)	Yes	UN-X approval	200 x 350 x 235	107 731
L	S 60 / 61	20 litres	PE-HD	natural (white)	Yes	UN-X approval	260 x 390 x 289	107 956
M	S 60 / 61	20 litres	PE-HD	natural (white)	No	with float	260 x 455 x 285	108 056
N	S 60 / 61	30 litres	PE-HD	natural (white)	Yes	UN-Y approval	290 x 400 x 380	107 959
O	S 60 / 61 & GL 45	30 litres	PE-HD	natural (white)	No	2 openings	290 x 400 x 380	108 115
P	S 70 / 71	60 litres	PE-HD	natural (white)	Yes	UN-Y approval	330 x 635 x 370	107 710
Q	S 90	10 litres	PE-HD	natural (white)	Yes	UN-Y approval	195 x 380 x 195	108 020
R	S 95	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	170 x 310 x 170	107 707
S	S 95	20 litres	PE-HD	natural (white)	No		260 x 390 x 290	107 733

SCAT Europe Overview

Range of Canisters



**THE No. 1 System
for Solvent Safety**



107 334



107 335



107 335



107 336

Product Group 3: Space-Saving Canisters

Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D	Part Number
A	S 50	5 litres	PP	natural (white)	No		65 x 335 x 330	107 998
B	S 50	5 litres	PP	natural (white)	No	with float	65 x 330 x 405	108 945

Product Group 4: Politainers

Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D	Part Number
C	GL 38	5 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 334	Box with UN-Y approval	178 x 178 x 178	107 330
D	GL 38	10 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 335	Box with UN-Y approval	228 x 228 x 228	107 331
E	S 60 / 61	10 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 335	Box with UN-Y approval	228 x 228 x 228	107 332
F	S 60 / 61	20 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 336	Box with UN-Y approval	285 x 285 x 285	107 333

SCAT Europe Overview

Range of Canisters



**THE No. 1 System
for Solvent Safety**



Product Group 5: Round Canisters

Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D	Part Number
A	B 53	2 litres	PE-HD	natural (white)	No		119 x 260 x 119	107 711
B	S 65	5 litres	PE-HD	natural (white)	No	with handle and measuring scale	167 x 330 x 167	107 722
C	S 65	10 litres	PE-HD	natural (white)	No	with handle and measuring scale	205 x 430 x 205	107 704
D	S 65	25 litres	PE-HD	natural (white)	No	w. carrying handles and measuring scale	278 x 580 x 278	107 720
-	S 65	60 litres	PE-HD	natural (white)	No	w. carrying handles and measuring scale	350 x 825 x 350	107 721
E	B 83	4 litres	PE-HD	natural (white)	No		155 x 338 x 155	107 712
F	B 83	10 litres	PP	natural (white)	No	with carrying handles	250 x 390 x 250	107 706
G	B 83	50 litres	PP	natural (white)	No	with carrying handles	380 x 680 x 380	107 730

Product Group 6: Flat Canisters

Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D	Part Number
H	S 70 / 71	20 litres	PE-HD	natural (white)	Yes	UN-Y approval	283 x 292 x 376	107 713

SCAT Europe Overview

Range of Laboratory Glass Bottles



**THE No. 1 System
for Solvent Safety**



Laboratory Glass Bottles

Image	Description	Thread	Capacity	Glass type	Form	Protective coating	Part Number
A	Original GL45 Screw Cap	GL 45	-	-	-	-	199 011
B	Laboratory bottle DURAN, GL45, 250 ml, Type 1	GL 45	250 ml	Clear glass	round	yes	101 994
C	Laboratory bottle DURAN, GL45, 250 ml, Type 3	GL 45	250 ml	Clear glass	round	-	501 117
D	Laboratory bottle DURAN, GL45, 250 ml, Type 2	GL 45	250 ml	Clear glass	square	-	501 112
E	Laboratory bottle DURAN, GL45, 500 ml, Type 1	GL 45	500 ml	Clear glass	round	yes	101 995
F	Laboratory bottle DURAN, GL45, 500 ml, Type 3	GL 45	500 ml	Clear glass	round	-	501 116
G	Laboratory bottle DURAN, GL45, 500 ml, Type 4	GL 45	500 ml	Brown glass	round	-	501 120
H	Laboratory bottle DURAN, GL45, 500 ml, Type 2	GL 45	500 ml	Clear glass	square	-	501 115
I	Laboratory bottle DURAN, GL45, 1 L, Type 1	GL 45	1 litre	Clear glass	round	yes	101 996
J	Laboratory bottle DURAN, GL45, 1 L, Type 2	GL 45	1 litre	Clear glass	round	-	501 113
K	Laboratory bottle DURAN, GL45, 1 L, Type 3	GL 45	1 litre	Brown glass	round	-	501 119
L	Laboratory bottle DURAN, GL45, 1 L	GL 45	1 litre	Clear glass	square	-	501 110

SCAT Europe Overview

Range of Laboratory Glass Bottles



**THE No. 1 System
for Solvent Safety**



Laboratory Glass Bottles

Image	Description	Thread	Capacity	Glass type	Form	Protective coating	Part Number
M	Laboratory bottle DURAN, GL45, 2 L, Type 1	GL 45	2 litres	Clear glass	round	yes	101 997
N	Laboratory bottle DURAN, GL45, 2 L, Type 2	GL 45	2 litres	Clear glass	round	-	501 118
O	Laboratory bottle DURAN, GL45, 5 L, Type 1	GL 45	5 litres	Clear glass	round	yes	101 998
-	Laboratory bottle DURAN, GL45, 5 L, Type 2	GL 45	5 litres	Clear glass	round	-	501 125
P	Laboratory bottle DURAN, GL45, 10 L	GL 45	10 litres	Clear glass	round	-	501 126
Q	Laboratory bottle DURAN, GLS80, 500 ml, Type 1	GLS80	500 ml	Clear glass	round	-	501 150
R	Laboratory bottle DURAN, GLS80, 1 L, Type 1	GLS80	1 litre	Clear glass	round	-	501 151
S	Laboratory bottle DURAN, GLS80, 1 L, Type 2	GLS80	1 litre	Brown glass	round	-	501 157
T	Laboratory bottle DURAN, GLS80, 2 L, Type 1	GLS80	2 litres	Clear glass	round	-	501 152



SCAT - Europe GmbH


Waldecker Straße 7
64546 Mörfelden-Walldorf
Germany

 **+49 (0) 61 05 / 30 55 86 - 0**

 **+49 (0) 61 05 / 30 55 86 - 99**

 **info@scat-europe.com**

 **www.scat-europe.com**



SCAT - Americas, Inc.

1704 East Blvd., Suite 101
Charlotte, NC 28203
USA

 **+1.704.438.9477**

 **+1.704.373.2603**

 **info@scat-americas.com**

 **www.scat-americas.com**

