

# Optimum Protection when collecting Liquid Waste.

In the hectic daily routine of a laboratory, things can easily go wrong. Open containers become a danger to people and the environment. Anyone who wants to get their waste disposal under control therefore needs a coherent concept. Our safety funnels are a one-time investment for many years of health and safety.

www.scat-europe.com





- Collect Liquid Waste Safely
- High-Quality Materials
- Robust and Durable
- Suitable for Organic Solvents

Level Co

Containers

Accessory

dendum

## Harmful Liquids Ever Spilled?

Open containers quickly create fire and health hazards. SCAT safety funnels can be safely closed and prevent ignition hazards. Also available in electrostatic conductive plastic!



# You know This? Comparison



## MARCO Safety Funnel with Hinged Lid

Smart protection when collecting liquid waste!



## **Safety Funnel MARCO**

#### **Technical Details**



#### **Safety Funnel MARCO**

#### **Ordering Information**























When our lid funnel MARCO is combined with the LISA, this combination goes by the name of JAN, see **page 94**.



Available in 150 mm and 190 mm length

### **Safety Funnel MARCO**

### **Ordering Information**





#### Informations about "delivery with adapter":

The safety funnel MARCO is available in various container screw-thread sizes. In order that we may produce in an environmentally-friendly and economic manner, certain variants are delivered together with a screw-thread adapter. Those to which this applies can be seen in the column "Delivery with Adapter".

Fig.	Part No.	Description	Screw-Thread of Container	Lance (220 mm)	Electrostatic conductive	Material	Level Control	Delivery with Adapter
A	318 952	Hinged lid Funnel MARCO	GL 45	-	-	PE-HD	-	-
	318 962	Hinged lid Funnel MARCO	GL 45	-	•	PE-HD-EL	-	-
B	318 992	Hinged lid Funnel MARCO	GL 45	•	-	PE-HD	-	-
	318 955	Hinged lid Funnel MARCO	S 50	-	-	PE-HD	-	•
	318 985	Hinged lid Funnel MARCO	S 50	•	•	PE-HD-EL	-	•
	318 995	Hinged lid Funnel MARCO	S 50	•	-	PE-HD	-	•
	318 953	Hinged lid Funnel MARCO	S 51	-	-	PE-HD	-	•
	318 983	Hinged lid Funnel MARCO	S 51	•	•	PE-HD-EL	-	•
	318 993	Hinged lid Funnel MARCO	S 51	•	-	PE-HD	-	•
	318 951	Hinged lid Funnel MARCO	S 55	-	-	PE-HD	-	-
	318 961	Hinged lid Funnel MARCO	S 55	-	•	PE-HD-EL	-	-
	318 981	Hinged lid Funnel MARCO	S 55	•	•	PE-HD-EL	-	-
	318 991	Hinged lid Funnel MARCO	S 55	•	-	PE-HD	-	-
	318 950	Hinged lid Funnel MARCO	S 60/61	-	-	PE-HD	-	-
G	318 960	Hinged lid Funnel MARCO	S 60/61	-	•	PE-HD-EL	-	-
D	318 980	Hinged lid Funnel MARCO	S 60/61	•	•	PE-HD-EL	-	-
	318 990	Hinged lid Funnel MARCO	S 60/61	•	-	PE-HD	-	-
8	317 637	Hinged lid Funnel MARCO	S 60/61	-	•	PE-HD-EL	• 150 mm	-
	317 651	Hinged lid Funnel MARCO	S 60/61	-	•	PE-HD-EL	• 190 mm	-
	318 954	Hinged lid Funnel MARCO	S 65	-	-	PE-HD	-	•
	318 964	Hinged lid Funnel MARCO	S 65	-	•	PE-HD-EL	-	•
	318 984	Hinged lid Funnel MARCO	S 65	•	•	PE-HD-EL	-	•
K	318 994	Hinged lid Funnel MARCO	S 65	•	-	PE-HD	-	•
<b>3</b>	318 999	Dirt sieve for hinged lid MARCO, PE-HD, blue						
G	318 998	Replacement lid for hinged lid funnel MARCO, PE-HD, blue						
	319 000	Dirt sieve for hinged lid MARCO, stainless steel, 105 mm, fine-mesh						
0	318 989	Dirt sieve for hinged lid MARCO, PE-HD-EL, black						
	318 988	Replacement lid for hinged lid funnel MARCO, PE-HD-EL, black						



**Disposal-Sets** 

Collect liquid HPLC waste safely: Sets including container and funnel, ready for immediate use - without time-consuming assembling of a configuration.





#### **Even more selection**

Further disposal sets, including electrostatic conductive versions, can be found in the "funnel sets" section of our website: **www.scat-europe.com** 

Scan the **QR code** to go directly to the funnel sets section of our website.



**Safety Funnel MARCO** 

### **Disposal-Sets**




Fig.	Part No.	Description	Contents		
A	307 316	Funnel Disposal-Set	Quantity	Description	See also
			1x	Hinged lid funnel MARCO, GL 45, PE-HD, with removable sieve	▶ Page 106
			1x	2,5 Litre canister, GL 45, PE-HD	▶ Page 136
	307 317	Funnel Disposal-Set	Quantity	Description	See also
			1x	Hinged lid funnel MARCO, GL 45, PE-HD, with removable sieve	▶ Page 106
			1x	5 Litre canister, GL 45, PE-HD	▶ Page 136
	307 318	Funnel Disposal-Set	Quantity	Description	See also
			1x	Hinged lid funnel MARCO, GL 45, PE-HD, with removable sieve	▶ Page 106
			1x	10 Litre canister, GL 45, PE-HD	▶ Page 136



#### **Safety Funnel ARNOLD**

#### **Technical Details**



# **Safety Funnel ARNOLD Ordering Information**













#### **Safety Funnel ARNOLD**

### **Ordering Information**





#### Informations about "delivery with adapter":

The safety funnel ARNOLD is available in various container screw-thread sizes. In order that we may produce in an environmentally-friendly and economic manner, certain variants are delivered together with a screw-thread adapter. Those to which this applies can be seen in the column "delivery with adapter".

Fig.	Part No.	Description	Screw-Thread of Container	Lance (220 mm)	Electrostatic conductive	Material	Delivery with Adapter
A	317 642	Ball-valve funnel ARNOLD	GL 45	•	-	PE-HD	-
В	317 622	Ball-valve funnel ARNOLD	GL 45	•	•	PE-HD-EL	-
	317 638	Ball-valve funnel ARNOLD	GL 45	-	•	PE-HD-EL	-
	317 649	Ball-valve funnel ARNOLD	S 50	•	-	PE-HD	•
	317 629	Ball-valve funnel ARNOLD	S 50	•	•	PE-HD-EL	•
	317 644	Ball-valve funnel ARNOLD	S 51	•	-	PE-HD	•
	317 624	Ball-valve funnel ARNOLD	S 51	•	•	PE-HD-EL	•
	317 632	Ball-valve funnel ARNOLD	B 53	•	-	PE-HD	•
	317 645	Ball-valve funnel ARNOLD	S 55	•	-	PE-HD	-
	317 625	Ball-valve funnel ARNOLD	S 55	•	•	PE-HD-EL	-
	317 641	Ball-valve funnel ARNOLD	S 60/61	•	-	PE-HD	-
	317 621	Ball-valve funnel ARNOLD	S 60/61	•	•	PE-HD-EL	-
	317 646	Ball-valve funnel ARNOLD	S 65	•	-	PE-HD	•
	317 626	Ball-valve funnel ARNOLD	S 65	•	•	PE-HD-EL	•
	317 648	Ball-valve funnel ARNOLD	S 70/71	•	-	PE-HD	•
	317 628	Ball-valve funnel ARNOLD	S 70/71	•	•	PE-HD-EL	•
	317 647	Ball-valve funnel ARNOLD	B 83	•	-	PE-HD	•
	317 627	Ball-valve funnel ARNOLD	B 83	•	•	PE-HD-EL	•
	317 643	Ball-valve funnel ARNOLD	S 90	•	-	PE-HD	•
	317 623	Ball-valve funnel ARNOLD	S 90	•	•	PE-HD-EL	•
G	317 635	Ball-valve funnel ARNOLD	S 95	•	•	PE-HD-EL	•
	317 630	Ball-valve funnel ARNOLD	Double thread R2" BSP/G2" + 2" TriSure	-	•	PE-HD-EL	-
D	319 001	Replacement lid for ball-valve funnel ARNOLD, PE-HD, blue					
<b>3</b>	317 640	Replacement sieve incl. splash protection for ball-valve funnel ARNOLD, PE-HD, blue					
<b>3</b>	319 002	Replacement lid for ball-valve funnel ARNOLD, PE-HD-EL, black					
G	317 620	Replacement sieve incl. splash protection for ball-valve funnel ARNOLD, PE-HD-EL, black					



**Disposal-Set** 

Collect liquid HPLC waste safely: Sets including container and funnel, ready for immediate use - without time-consuming assembling of a configuration.





#### **Even more selection**

Further disposal sets, including electrostatic conductive versions, can be found in the "funnel sets" section of our website: **www.scat-europe.com** 

Scan the **QR code** to go directly to the funnel sets section of our website.



# **Safety Funnel ARNOLD Disposal-Set**



Fig.	Part No.	Description	Contents		
A	<b>307 450</b> Funnel Disposal-Set		Quantity	Description	See also
			1x	Safety Funnel with ball-valve ARNOLD, S 55, PE-HD, sieve & splash protection	▶ Page 113
			1x	10 Litre canister S 55, PE-HD	▶ Page 137

# Funnel ARNOLD for Barrels

Double thread R2" BSP/G2" (m) + BCS 56x4 (m)

#### **Version without Ball-Valve**

- ARNOLD funnels without ball valve are adapted for high flow rates or highly viscous liquids
- Equipped with an additional, removable dirt sieve (instead of ball-valve)





#### **Double Thread**

R2" BSP/G2" (m) + BCS 56x4 (m)

Fig.	Part No.	Description	Material
A	317 631	Barrel-Funnel ARNOLD with lid, double thread R2" BSP/G2" (m) + BCS 56x4 (m), with splash protection, removable sieve and grounding cable, 2 sieves	PE-HD-EL
	317 630	Barrel-Funnel ARNOLD <b>with ball-valve</b> and lid, double thread R2" BSP/G2" (m) + BCS 56x4 (m), with splash protection, removable sieve and grounding cable	PE-HD-EL

# Funnel ARNOLD with Level Control

For safe filling of barrels and large containers









Fig.	Part No.	Description	Material
A	317 636	Funnel ARNOLD GL 45 with lid, 2 sieves	PE-HD
В	317 633	Funnel ARNOLD GL 45 with lid, 2 sieves	PE-HD-EL
G	107 892	Adapter for funnels with mechanical level control, barrel double thread R2" BSP/G2" (m) to GL 45 (m)	PP
D	107 885 Adapter for funnels with mechanical level control, barrel double thread R2" BSP/G2" (m) to GL 45 (m)		PE-HD-EL
	317 660	Set consisting of <b>A</b> and <b>G</b>	-
	317 650	Set consisting of <b>B</b> and <b>D</b>	-

