



Maximize Your Workflow Efficiency and Your Budget

Agilent ValueLab Filter and Collection Plates

Exceptional Versatility

Our comprehensive selection of plates is designed to meet a broad spectrum of applications, allowing you to stay flexible in your processes.

Reliable Workflow Optimization

ValueLab filter and collection plates help streamline your sample preparation, ensuring efficiency in your lab operations.

Performance Meets Affordability

Agilent ValueLab products are tailored to meet your needs, delivering the right product for the right application—at the right price.

Choose From a Variety of Cost-Effective Solutions

Discover the exceptional value and performance of Agilent's ValueLab filter and collection plates. Engineered to meet the demands of a diverse range of applications, our products offer unparalleled versatility without compromising your budget. With ValueLab, you can trust that you're investing in a reliable solution designed to optimize your workflow efficiency. Whether you're working in research, diagnostics, or any other field, our filter and collection plates ensure consistent, high-quality results.

Ordering Information

Full Skirt Filter Plates

Wells	Filter	Pore Size	Volume	Quantity	Part Number
96	PE-frit	20 µm	1 mL	50/pk	5610-2204
96	PE-frit	5 µm	1 mL	50/pk	5610-2205
96	CA	0.45 µm	0.35 mL	50/pk	5610-2206
96	PVDF	0.45 µm	0.35 mL	50/pk	5610-2207
96	PVDF	0.22 µm	0.35 mL	50/pk	5610-2208
96	PTFE	0.45 µm	0.35 mL	50/pk	5610-2209
96	PTFE	0.22 µm	0.35 mL	50/pk	5610-2210
96	PTFE-HI	0.45 µm	0.35 mL	50/pk	5610-2211
96	PTFE-HI	0.22 µm	0.35 mL	50/pk	5610-2212
96	PTFE*	0.22 µm	1.8 mL	50/pk	5610-2213

* Full skirt protein precipitation plate



Collection Plates

Wells	Well Shape	Volume	Quantity	Part Number
96	Round	0.35 mL	100/pk	5610-2196
96	Round	0.5 mL	100/pk	5610-2197
96	Round	1.1 mL	100/pk	5610-2198
96	Square	2.3 mL	100/pk	5610-2199
48	Square	4.6 mL	100/pk	5610-2202

Sealing Gaskets

Wells	Well Shape	Quantity	Part Number
96	Round	100/pk	5610-2200
96	Square	100/pk	5610-2201
48	Square	100/pk	5610-2203



Related Products for Sample Filtration

Captiva Well Plates, Filter Plates, and Cartridges

Captiva's portfolio of sample filtration solutions is dedicated to extending the lifespan of your columns, reducing instrument downtime, and ensuring sample purity. These filtration products effectively eliminate precipitated proteins and offer robust safeguards against column clogging, resulting in the seamless and efficient operation of your laboratory.

Learn more:

www.agilent.com/chem/captiva

DE58490788

This information is subject to change without notice.

© Agilent Technologies, Inc. 2023
Published in the USA, November 1, 2023
5994-6823EN