


More Choices, Less Stress


Agilent InfinityLab well plates
and sealing mats



Choosing the right well plates can help you avoid such pitfalls as missed data points, sample rerun, and a failure to meet project deadlines.

New Agilent InfinityLab well plates for the HPLC workflow can help you produce the highest-quality data for your most common screening challenges. Our comprehensive portfolio:

- Features the right formats to ensure that your autosampler setup is perfectly equipped to deal with diverse HPLC applications
- Complies with worldwide ANSI and SLAS standards
- Offers different packaging quantities to suit your lab size and needs
- Perfect fit into your InfinityLab workflow, so you can choose from a comprehensive selection of well plates to complete your analysis



Pressed for time? The Agilent online store is a fast, easy way to order—and reorder—your well plates.

At the Agilent online store, well plates and other consumables are just a few clicks away. What's more, our "favorites" tab makes reordering as easy as 1-2-3:

1. Click "My favorites"
2. Click "Create new list"
3. Add part numbers

You can even create separate lists for different workflows, instruments, or applications.

Try it now: www.agilent.com/chem/store

Ordering Information

Well plates

All well plates are made in polypropylene

No. of wells	Well volume (mL)	Well shape	Bottom shape	Height	Units per pack	Well plate part no.	Recommended mat
96	2.2	Square	U	41 mm	30	5043-9300	5043-9319
96	2.0	Round	U	45 mm	30	5043-9302	5043-9317 / 5043-9318
96	1.0	Round	U	32 mm	50	5043-9305	5043-9317 / 5043-9318
96	1.2	Round	U	27 mm	25	5043-9308	5043-9317 / 5043-9318
96	1.2	Round	U	27 mm	50	5043-9309	5043-9317 / 5043-9318
96	0.500	Round	U	14 mm	30	5043-9310	5043-9317 / 5043-9318
96	0.500	Round	U	14 mm	120	5043-9311	5043-9317 / 5043-9318
96	0.330	Round	V	14 mm	25	5043-9312	5043-9317 / 5043-9318
96	0.330	Round	V	14 mm	50	5043-9313	5043-9317 / 5043-9318
96	0.330	Round	V	14 mm	100	5043-9314	5043-9317 / 5043-9318
384	0.190	Square	V	22 mm	25	5043-9315	5043-9320

Mats

All mats are made in silicone and are pre-slitted for best usage with HPLC samplers

Part no.	No. of wells	Well shape	Units per pack
5043-9317	96	Round	50
5043-9318	96	Round	100
5043-9319	96	Square	50
5043-9320	384	Square	50



Agilent InfinityLab is an optimized portfolio of LC and LC/MS instruments, columns, and supplies designed to work together in perfect harmony.

Combined with Agilent OpenLab software and Agilent CrossLab services, it provides you with end-to-end solutions and support to make every day more productive.

Find the right well plate for your application:

www.agilent.com/chem/well-plates

Place your order quickly and easily:

www.agilent.com/chem/store

For Research Use Only. Not for use in diagnostic procedures.

This information is subject to change without notice.

© Agilent Technologies, Inc. 2018
Published in the USA, November 1, 2018
5994-0145EN